AVANOS\*


# MIC\* GASTROSTOMY FEEDING TUBE


WITH ENFit® CONNECTORS

# MIC\* BOLUS GASTROSTOMY FEEDING TUBE

WITH ENFit® CONNECTOR

Instructions for Use

Diameter	Balloon Volume	Not made with natural rubber latex	Product is NOT made with DEHP as a plasticizer
Single Use Only	<b>STERILE</b> Sterilized by Gamma Irradiation	Do not resterilize	Do not use if package is damaged
Rx Only	<b>MR</b> MR Safe	Caution	Consult instructions for use

## Instructions for Use

Rx Only: Federal Law (USA) restricts this device to sale by or on the order of a physician.

## Description

The AVANOS® family of MIC® gastrostomy feeding tubes (Fig. 1 & 2) allows for delivery of enteral nutrition and medication directly into the stomach and/or gastric decompression.

## Indications for Use

The AVANOS® family of MIC® gastrostomy feeding tubes is indicated for use in patients who require long term feeding, are unable to tolerate oral feeding, who are at low risk for aspiration, require gastric decompression and/or medication delivery directly into the stomach.

## Contraindications

Contraindications for placement of a gastrostomy feeding tube include, but are not limited to:

- Colonic interposition
- Ascites
- Portal hypertension
- Peritonitis
- Uncorrected coagulopathy
- Uncertainty as to gastrostomy tract direction and length (abdominal wall thickness)
- Lack of adherence of the stomach to the abdominal wall (replacement only)
- Lack of established gastrostomy tract (replacement only)
- Evidence of infection around stoma site (replacement only)
- Presence of multiple stoma fistulous tracts (replacement only)

## ⚠ Warning

**Do not reuse, reprocess, or sterilize this medical device. Reuse, reprocessing, or sterilization may 1) adversely affect the known biocompatibility characteristics of the device, 2) compromise the structural integrity of the device, 3) lead to the device not performing as intended, or 4) create a risk of contamination and cause the transmission of infectious diseases resulting in patient injury, illness, or death.**

## Complications

The following complications may be associated with any gastrostomy feeding tube:

- Nausea, vomiting, abdominal bloating or diarrhea
- Aspiration
- Peristomal pain
- Abscess, wound infection and skin breakdown
- Pressure necrosis
- Hypergranulation tissue
- Intraperitoneal leakage
- Buried bumper syndrome
- Peristomal leakage
- Balloon failure or tube dislodgement
- Tube clog
- Gastrointestinal bleeding and/or ulceration
- Gastric outlet obstruction
- Ileus or gastroparesis
- Bowel and gastric volvulus

Other complications such as abdominal organ injury may be associated with the procedure to place the feeding tube.

## Placement

The AVANOS® family of gastrostomy feeding tubes may be placed in the following ways:

- Ⓢ Surgically
- Ⓡ Percutaneously under fluoroscopic (radiologic) guidance
- ⓔ Percutaneously under endoscopic guidance
- ⓧ A replacement to an existing device using an established stoma tract

**⚠ Warning: A gastrostomy must be performed to affix the stomach to the anterior abdominal wall, the feeding tube insertion site identified and stoma tract dilated prior to initial tube insertion to ensure patient safety and comfort.**

**⚠ Warning: Do not use the retention balloon of the feeding tube as a gastrostomy device. The balloon may burst and fail to attach the stomach to the anterior abdominal wall.**

**⚠ Warning: The insertion site for infants and children should be high on the greater curvature to prevent occlusion of the pylorus when the balloon is inflated.**

## Ⓢ Ⓡ ⓔ ⓧ Tube Preparation

**⚠ Warning: Verify package integrity. Do not use if package is damaged or sterile barrier compromised.**

1. Select the appropriate size MIC® gastrostomy feeding tube, remove from the package and inspect for damage.
2. Using a male Luer syringe, inflate the balloon (Fig. 1-E & 2-D) with water through the Balloon Inflation Port (Fig. 1-A & 2-A). Do not use air.
3. Inflate the balloon with 2-3 ml of water for 12 Fr low volume tubes identified by LV following the REF code number.
4. Inflate the balloon with 3-5 ml of water for 14 Fr and 16 Fr low volume tubes identified by LV following the REF code number.
5. Inflate the balloon with 7-10 ml of water for standard tubes.
6. Remove the syringe and verify balloon integrity by gently squeezing the balloon to check for leaks. Visually inspect the balloon to verify symmetry. Symmetry may be achieved by gently rolling the balloon between the fingers. Reinsert the syringe and remove all the water from the balloon.

4. Check the external retention bolster (Fig. 1-D & 2-C). The bolster should slide along the tube with moderate resistance.
5. Inspect the entire length of the tube for any irregularities.
6. Lubricate the tip of the tube with a water soluble lubricant. Do not use mineral oil. Do not use petroleum jelly.

## Ⓢ Suggested Surgical Procedure (Stamm Gastrostomy)

1. Through a midline laparotomy, identify the pylorus and the superior epigastric artery in the abdominal wall.
2. The gastrostomy site should be 10-15 cm from the pylorus on the greater curvature of the stomach. The gastrostomy site should also be at least 3 cm away from the costal margin to prevent damage to the retention balloon by abrasions during movement.
3. Place two concentric purse string sutures around the site. Leave the purse string needles in place.
4. On the anterior parietal peritoneum, select an exit site that approximates the gastrostomy. Avoid the superior epigastric artery, drains, or other stomas.
5. Make a stab wound from the anterior parietal peritoneum to the extra abdominal surface. Insert the tube from the outside to the inside the abdominal cavity.
 

**Note:** The use of a right angle clamp may be used to facilitate placement.

**⚠ Caution:** Never use a clamp with sharp teeth or tenaculum to pull the tube into position. This will damage the tube.
6. Using two Babcock clamps on the anterior stomach surface, "tent" the stomach.
7. Use electrocautery or a scalpel to open the stomach.
8. Dilate the enterotomy with a hemostat.

## Ⓡ ⓔ Suggested Site Preparation

1. Use standard Radiologic or Endoscopic techniques to visualize and prepare for gastrostomy tube placement.
2. Verify that no abnormalities are present that could pose a contraindication to placement of the tube and place the patient in the supine position.
3. Select a gastrostomy site that is free of major vessels, viscera, and scar tissue. The site is usually one third the distance from the umbilicus to the left costal margin at the midclavicular line.
4. Prep and drape the selected insertion site per facility protocol.

## Ⓡ ⓔ Gastrostomy Placement

**⚠ Warning: A three point gastrostomy in a triangle configuration is recommended to ensure attachment of the gastric wall to the anterior abdominal wall.**

1. Place a skin mark at the tube insertion site. Define the gastrostomy pattern by placing three skin marks equidistant from the tube insertion site and in a triangle configuration.
 

**⚠ Caution:** Allow adequate distance between the insertion site and gastrostomy placement to prevent interference of the T-Fastener and inflated balloon.
2. Localize the puncture sites with 1% lidocaine and administer local anesthesia to the skin and peritoneum.
3. Place the first T-Fastener and confirm intragastric position. Repeat the procedure until all three T-Fasteners are inserted at the corners of the triangle.
4. Secure the stomach to the anterior abdominal wall and complete the procedure.

## Ⓡ ⓔ Create the Stoma Tract

1. Create the stoma tract with the stomach still insufflated and in apposition to the abdominal wall. Identify the puncture site at the center of the gastrostomy pattern. With fluoroscopic or endoscopic guidance, confirm that the site overlies the distal body of the stomach below the costal margin and above the transverse colon.

**⚠ Warning: Avoid the epigastric artery that courses at the junction of the medial two-thirds and lateral one-third of the rectus muscle.**

2. Anesthetize the puncture site with local injection of 1% lidocaine down to the peritoneal surface.
3. Insert a .038" compatible introducer needle at the center of the gastrostomy pattern into the gastric lumen.

**⚠ Warning: Take care not to advance the puncture needle too deeply in order to avoid puncturing the posterior gastric wall, pancreas, left kidney, aorta or spleen.**

**Note:** For gastrostomy tube placement, the best angle of insertion is a true right angle to the surface of the skin. The needle, however, should be directed toward the pylorus if conversion to a jejunal feeding tube is anticipated.

4. Use fluoroscopic or endoscopic visualization to verify correct needle placement. Additionally, to aid in verification, a water filled syringe may be attached to the needle hub and air aspirated from the gastric lumen.
 

**Ⓡ Note:** Contrast may be injected upon return of air to visualize gastric folds and confirm position.
5. Advance a J tip guidewire, up to .038", through the needle and into stomach. Confirm position.
6. Remove the introducer needle, keeping the J tip guidewire in place and dispose of according to facility protocol.

## Ⓡ ⓔ Dilation

1. Use a #11 scalpel blade to create a skin incision that extends alongside the guidewire, downward through the subcutaneous tissue and fascia of the abdominal musculature. After the incision is made, dispose of the scalpel according to facility protocol.
2. Advance a dilator over the guidewire and dilate the stoma tract at least four French sizes larger than the enteral feeding tube being placed.
3. Remove the dilator over the guidewire, leaving the guidewire in place.

**Note:** After dilation, a peel-away sheath may be used to facilitate advancement of the tube through the stoma tract.

## Ⓢ Tube Placement

- Advance the MIC® gastrostomy feeding tube until the balloon is in the stomach.
- Using a male Luer syringe, inflate the balloon.
  - Inflate 12 Fr LV balloons with 2-3 ml of sterile or distilled water.
  - Inflate 14 Fr and 16 Fr LV balloons with 3-5 ml of sterile or distilled water.
  - Inflate the standard balloon with 7-10 ml of sterile or distilled water.

**⚠️ Caution:** Do not exceed 5 ml total balloon volume in the 12 Fr LV balloon, 7 ml total balloon volume in the 14 Fr and 16 Fr LV balloons, and 15 ml total balloon volume in the Standard balloon. Do not use air. Do not inject contrast into the balloon.
- Tie the purse string sutures around the tube.
- Gently pull the tube up and away from the abdomen until the balloon contacts the inner stomach wall.
- Use the purse string sutures to attach the stomach to the peritoneum. Take care to avoid puncture of the balloon.
- Clean the residual fluid or lubricant from the tube and stoma.
- Gently slide the external retention bolster to approximately 2-3 mm above the skin. Do not suture the bolster to the skin.

## Ⓡⓔ Tube Placement

- Advance the distal end of the tube over the guidewire, through the stoma tract and into the stomach.
- Verify that the tube is in the stomach, remove the guidewire or peel-away sheath if utilized and inflate the balloon.
- Using a male Luer syringe, inflate the balloon.
  - Inflate 12 Fr LV balloons with 2-3 ml of sterile or distilled water.
  - Inflate 14 Fr and 16 Fr LV balloons with 3-5 ml of sterile or distilled water.
  - Inflate the standard balloon with 7-10 ml of sterile or distilled water.

**⚠️ Caution:** Do not exceed 5 ml total balloon volume in the 12 Fr LV balloon, 7 ml total balloon volume in the 14 Fr and 16 Fr LV balloons, and 15 ml total balloon volume in the Standard balloon. Do not use air. Do not inject contrast into the balloon.
- Gently pull the tube up and away from the abdomen until the balloon contacts the inner stomach wall.
- Clean the residual fluid or lubricant from the tube and stoma.
- Gently slide the external retention bolster to approximately 2-3 mm above the skin. Do not suture the bolster to the skin.

## ⓈⓇⓔⓧ Verify Tube Position and Patency

- Attach an ENFit® syringe with 10 ml water to an access port (Fig. 1-B, 1-C, & 2-B) of the gastrostomy feeding tube. Aspirate gastric contents. When air or gastric contents are observed, flush the tube.
- Check for moisture around the stoma. If there are signs of gastric leakage, check the tube position and placement of the external retention bolster. Add sterile or distilled water as needed in 1-2 ml increments. Do not exceed balloon capacity as indicated previously.
- Begin feeding only after confirmation of proper patency, placement and according to physician instructions.

## ⓧ Tube Removal

- Ensure that this type of tube can be replaced at the bedside.
- Assemble all equipment and supplies, cleanse hands using aseptic technique and apply clean, powder-free gloves.
- Rotate the tube 360 degrees to ensure the tube moves freely and easily.
- Firmly insert a male Luer syringe into the balloon port and withdraw all the fluid from the balloon.
- Apply counter pressure to the abdomen and remove the tube with gentle, but firm traction.

**Note:** If resistance is encountered, lubricate the tube and stoma with water soluble lubricant. Simultaneously push and rotate the tube. Gently manipulate the tube free. If the tube will not come out, refill the balloon with the prescribed amount of water and notify the physician. Never use excessive force to remove a tube.

**⚠️ Caution:** Never attempt to change the tube unless trained by the physician or other health care provider.

**⚠️ Warning:** After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

## ⓧ Replacement Procedure

- Cleanse the skin around the stoma site and allow the area to air dry.
- Select the appropriate size Gastrostomy feeding tube and prepare according to the instructions in the Tube Preparation section above.
- Gently insert the Gastrostomy through the stoma into the stomach.
- Using a male Luer syringe, inflate the balloon.
  - Inflate 12 Fr LV balloons with 2-3 ml of sterile or distilled water.
  - Inflate 14 Fr and 16 Fr LV balloons with 3-5 ml of sterile or distilled water.
  - Inflate the Standard balloon with 7-10 ml of sterile or distilled water.

**⚠️ Caution:** Do not exceed 5 ml total balloon volume in the 12 Fr LV balloon, 7 ml total balloon volume in the 14 Fr and 16 Fr LV balloons, and 15 ml balloon volume in the Standard balloon. Do not use air. Do not inject contrast into the balloon.
- Gently pull the tube up and away from the abdomen until the balloon contacts the inner stomach wall.
- Clean the residual fluid or lubricant from the tube and stoma.
- Gently slide the external retention bolster to approximately 1-2 mm above the skin.
- Verify proper tube position according to the instructions in the Verify Tube Position section above.

## Tube Patency Guidelines

Proper tube flushing is the best way to avoid clogging and maintain tube patency. The following are guidelines to avoid clogging and maintain tube patency.

- Flush the feeding tube with water every 4-6 hours during continuous feeding, anytime the feeding is interrupted, before and after every intermittent feeding, or at least every 8 hours if the tube is not being used.
- Flush the feeding tube after checking gastric residuals.
- Flush the feeding tube before and after medication administration and between medications. This will prevent the medication from interacting with formula and potentially causing the tube to clog.
- Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and to mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in warm water before administering through the feeding tube. Never crush enteric-coated medication or mix medication with formula.
- Avoid using acidic fluids such as cranberry juice and cola beverages to flush feeding tubes as the acidic quality when combined with formula proteins may actually contribute to tube clogging.

## General Flushing Guidelines

Flush the feeding tube with water using an ENFit® syringe every 4-6 hours during continuous feeding, anytime the feeding is interrupted, at least every 8 hours if the tube is not being used, or per clinician's instructions. Flush the feeding tube after checking gastric residuals. Flush the feeding tube before and after each medication administration. Avoid using acidic irrigants such as cranberry juice and cola beverages to flush feeding tubes.

- Use a 30 to 60 ml ENFit® syringe. Do not use smaller size syringes as this can increase pressure on the tube and potentially rupture smaller tubes.
- Ensure the second access port (if applicable) is closed with the tethered cap prior to flushing.
- Use room temperature water for tube flushing. Sterile water may be appropriate where the quality of municipal water supplies is of concern. The amount of water will depend on the patient's needs, clinical condition, and type of tube, but the average volume ranges from 10 to 50 ml for adults, and 3 to 10 ml for infants. Hydration status also influences the volume used for flushing feeding tubes. In many cases, increasing the flushing volume can avoid the need for supplemental intravenous fluid. However, individuals with renal failure and other fluid restrictions should receive the minimum flushing volume necessary to maintain patency.
- Do not use excessive force to flush the tube. Excessive force can perforate the tube and can cause injury to the gastrointestinal tract.
- Document the time and amount of water used in the patient's record. This will enable all caregivers to monitor the patient's needs more accurately.

## Nutrition Administration

- Open the cap to an access port of the gastrostomy tube.
- Use an ENFit® syringe to flush the tube with the prescribed amount of water as described in the General Flushing Guidelines.
- Remove the flushing syringe from the access port.
- Securely connect an ENFit® feed set or an ENFit® syringe to the access port.

**⚠️ Caution:** Do not over-tighten the feed set connector or the syringe to the access port.
- Complete feeding per the clinician's instructions.
- Remove the feed set or syringe from the access port.
- Use an ENFit® syringe to flush the tube with the prescribed amount of water as described in the General Flushing Guidelines.
- Remove the flushing syringe from the access port.
- Close the cap to the access port.

## Medication Administration

Use liquid medication when possible and consult the pharmacist to determine if it is safe to crush solid medication and mix with water. If safe, pulverize the solid medication into a fine powder form and dissolve the powder in warm water before administering through the feeding tube. Never crush enteric coated medication or mix medication with formula.

- Open the cap to an access port of the gastrostomy tube.
- Use an ENFit® syringe to flush the tube with the prescribed amount of water as described in the General Flushing Guidelines.
- Remove the flushing syringe from the access port.
- Securely connect an ENFit® syringe containing the medication to the access port.

**⚠️ Caution:** Do not over-tighten the syringe to the access port.
- Deliver the medication by depressing the ENFit® syringe plunger.
- Remove the syringe from the access port.
- Use an ENFit® syringe to flush the tube with the prescribed amount of water as described in the General Flushing Guidelines.
- Remove the flushing syringe from the access port.
- Close the cap to the access port.

## Gastric Decompression

Gastric decompression may be performed via either gravity drainage or low intermittent suction.

- Open the cap to an access port of the gastrostomy tube.
- For gravity drainage, place the opened access port of the gastrostomy tube directly over the opening of an appropriate container.

**Note:** Ensure the open access port is positioned below the stoma.
- For low intermittent suction, connect an ENFit® syringe to the access port.
- Apply low intermittent suction by slowly retracting the plunger of the syringe in short intervals.

**⚠️ Warning:** Do not use continuous or high pressure suction. High pressure could collapse the tube or injure the stomach tissue and cause bleeding.
- Disconnect the decompression syringe from the access port.
- Use an ENFit® syringe to flush the tube with the prescribed amount of water as described in the General Flushing Guidelines.

7. Remove the flushing syringe from the access port.
8. Close the cap to the access port.

## Balloon Maintenance

Precise balloon life cannot be predicted. Silicone balloons generally last 1–8 months, but the life span of the balloon varies according to several factors. These factors may include medications, volume of water used to inflate the balloon, gastric pH, and tube care.

Check the water volume in the balloon once a week.

- Insert a male Luer syringe into the Balloon Inflation Port (BAL.) and withdraw the fluid while holding the tube in place. Compare the amount of water in the syringe to the amount recommended or the amount initially prescribed and documented in the patient record. If the amount is less than recommended or prescribed, refill the balloon with the water initially removed, then draw up and add the amount needed to bring the balloon volume up to the recommended and prescribed amount of water. Be aware as you deflate the balloon there may be some gastric contents that can leak from around the tube. Document the fluid volume, the amount of volume to be replaced (if any), the date and time.
- Wait 10–20 minutes and repeat the procedure. The balloon is leaking if it has lost fluid, and the tube should be replaced. A deflated or ruptured balloon could cause the tube to dislodge or be displaced. If the balloon is ruptured, it will need to be replaced. Secure the tube into position using tape, then follow facility protocol and/or call the physician for instructions.  
**⚠️ Caution:** Refill the balloon using sterile or distilled water, not air or saline. Saline can crystallize and clog the balloon valve or lumen, and air may seep out and cause the balloon to collapse. Be sure to use the recommended amount of water as over-inflation can obstruct the lumen or decrease balloon life and under-inflation will not secure the tube properly.

## Daily Care & Maintenance Check List

- Assess the patient for any signs of pain, pressure or discomfort.
- Assess the stoma site for any signs of infection, such as redness, irritation, edema, swelling, tenderness, warmth, rashes, purulent, or gastrointestinal drainage. Assess the patient for any signs of pressure necrosis, skin breakdown, or hypergranulation tissue.
- Clean the stoma site using warm water and mild soap.
- Use a circular motion moving from the tube outwards.
- Rinse thoroughly and dry well.
- Assess the tube for any abnormalities such as damage, clogging, or abnormal discoloration.
- Clean the feeding tube using warm water and mild soap being careful not to pull or manipulate the tube excessively.
- Rinse thoroughly and dry well.
- Clean the gastric and balloon inflation ports. Use a cotton tip applicator or soft cloth to remove all residual formula and medication.
- Rotate the tube 360 degrees plus a quarter turn daily.
- Verify that the external bolster rests 1–2 mm above the skin.
- Flush the feeding tube as described in the General Flushing Guidelines section above.

## Tube Occlusion

Tube occlusion is generally caused by:

- Poor flushing techniques
- Failure to flush after measurement of gastric residuals
- Inappropriate administration of medication
- Pill fragments
- Thick formulas, such as concentrated, blenderized, or enriched formulas that are generally thicker and may contain particulates
- Formula contamination that leads to coagulation
- Reflux of gastric or intestinal contents up the tube

## To Unclog a Tube

1. Make sure that the feeding tube is not kinked or clamped off.
2. If the clog is visible above the skin surface, gently massage or milk the tube between fingers to break up the clog.
3. Connect a 30 to 60 ml ENFit® syringe filled with warm water into the appropriate access port of the tube and gently pull back on then depress the plunger to dislodge the clog.
4. If the clog remains, repeat step #3. Gentle suction alternating with syringe pressure will relieve most obstructions.
5. If this fails, consult with the physician. Do not use cranberry juice, cola drinks, meat tenderizer or chymotrypsin, as they can actually cause clogs or create adverse reactions in some patients. If the clog is stubborn and cannot be removed, the tube will have to be replaced.

**⚠️ Caution:** Do not insert foreign objects through the tube.


## MRI Safety Information

The MIC® Gastrostomy Feeding Tubes are MR Safe.

**⚠️ Warning:** For enteral nutrition and/or enteral medication only.

For more information, please call 1-844-4AVANOS (1-844-4AVANOS (1-844-428-2667)) in the United States, or visit our web site at [www.avanos.com](http://www.avanos.com).

Educational Booklets: "A Guide to Proper Care" and "A Stoma Site and Enteral Feeding Tube Troubleshooting Guide" is available upon request. Please contact your local representative or contact Customer Care.

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Distributed in the USA by Avanos Medical Sales, LLC,  
5405 Windward Parkway, Alpharetta, GA 30004 USA.  
In USA, 1-844-4AVANOS (1-844-428-2667). [www.avanos.com](http://www.avanos.com)  
Avanos Medical, Inc., 5405 Windward Parkway, Alpharetta, GA 30004 USA.  
 Avanos Medical Belgium BVBA, Leonardo Da Vincilaan 1,  
1930 Zaventem, Belgium.  
Sponsored in Australia by Avanos Medical Australia Pty Ltd,  
475 Victoria Avenue, Chatswood, NSW 2067, Australia.

\*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.  
© 2018 AVNS. All rights reserved. 2018-12-20  
15-M1-276-01 / 70210674