
# Shiley™

## Tracheostomy Tube

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
it	Istruzioni per l'uso
es	Instrucciones de uso
pt	Instruções de uso
nl	Gebruiksaanwijzing
sv	Bruksanvisning
da	Brugsanvisning
no	Bruksanvisning
fi	Käyttöohjeet
el	Οδηγίες χρήσης
zh	使用说明

1. LPC & FEN

					
Product Size	I.D. (mm)	O.D. (mm)	Length (mm)	Leak Test Inflation Volume (cc)	Cuff Resting Diameter (mm)
4	5.0	9.4	65.0	9	20
6	6.4	10.8	76.0	12	24
8	7.6	12.2	81.0	15	27
10	8.9	13.8	81.0	18	29







2. CFS & CFN

			
Product Size	I.D. (mm)	O.D. (mm)	Length (mm)
4	5.0	9.4	65.0
6	6.4	10.8	76.0
8	7.6	12.2	81.0
10	8.9	13.8	81.0

3. LGT

			
Product Size	I.D. (mm)	O.D. (mm)	Length (mm)
6	6.4	11.1	50.0
8	7.6	12.6	50.0
10	8.9	13.7	50.0

4.

	LPC	FEN	CFS	CFN	LGT
	X	X	X	X	X
	X	X	X	X	X
		X		X	
		X		X	
			X	X	X(2)
			X		
	X	X	X	X	X

Identification of a substance that is contained or present within the product or packaging.

Identification of a substance that is not contained or present within the product or packaging.

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This product contains DEHP. When used as indicated, very limited exposure to trace amounts of DEHP may occur. There is no clear clinical evidence that this degree of exposure increases clinical risk. However, in order to minimize risk of DEHP exposure in children and nursing or pregnant women, this product should only be used as directed.

**This instruction insert applies to the following Shiley™ Tracheostomy products:**

LPC	Tracheostomy Tube Cuffed with Inner Cannula
FEN	Tracheostomy Tube Cuffed with Inner Cannula Fenestrated
CFS	Tracheostomy Tube Cuffless with Inner Cannula
CFN	Tracheostomy Tube Cuffless with Inner Cannula Fenestrated
LGT	Laryngectomy Tube

#### Instructions for Use

##### **WARNINGS:**

**CAREFULLY READ THESE INSTRUCTIONS PRIOR TO USE. THESE INSTRUCTIONS ARE APPLICABLE ONLY TO THE SHILEY DOUBLE CANNULA TRACHEOSTOMY TUBES (LPC, FEN, CFS, CFN) AND SHILEY LARYNGECTOMY TUBE (LGT) WITH REUSEABLE INNER CANNULAE.**

**STERILE ONLY IF THE PROTECTIVE BAG IS NOT OPEN, DAMAGED OR BROKEN. DO NOT RESTERILIZE. DO NOT EXPOSE TO TEMPERATURES ABOVE 120°F (49°C).**

**WHEN USING A LASER IN CONJUNCTION WITH THESE DEVICES, CARE MUST BE TAKEN TO AVOID CONTACT OF THE LASER BEAM WITH THE TUBE. SUCH CONTACT, ESPECIALLY IN THE PRESENCE OF OXYGEN-ENRICHED MIXTURES, COULD RESULT IN RAPID COMBUSTION OF THE TUBE WITH HARMFUL THERMAL EFFECTS AND WITH EMISSION OF CORROSIVE AND TOXIC COMBUSTION PRODUCTS, INCLUDING HYDROCHLORIC ACID (HCl).**

**DURING AN MRI SCAN THE PILOT BALLOON SHOULD BE SECURED NEAR THE Y CONNECTOR OF THE VENTILATOR CIRCUIT AT LEAST 3 CM FROM THE AREA OF INTEREST TO PREVENT MOVEMENT AND IMAGE DISTORTION.**

## POST FOR USE BY ALL TRACHEOSTOMY CARE PERSONNEL

**CAUTION:** Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

**CAUTIONS:**

*The Shiley™ Tracheostomy Tube and Accessory Products are composed of medical grade plastics which do not contain latex material.*

*Precautions should be taken when discarding any of the tracheostomy tube devices and accessories. Disposal of medical devices should be made in accordance with applicable national regulations for biologically hazardous waste.*

Contents of Packages - Refer to Table 4.

Description - Refer to Tables 1 through 3 For Sizing Dimensions

**NOTE:** The sizing dimensions as noted in Tables 1 through 3, refer to the I.D. as the inside diameter of the inner cannula at the smallest location. O.D. is the outside diameter of the outer cannula (not including cuff features). Length is the distance from the neck plate to the distal tip on the tube's centerline.

Shiley Tracheostomy Tubes are double cannula tracheostomy tubes with reusable inner cannulae (LPC, FEN, CFS, CFN, LGT) and twist-lock connectors. These Shiley Tracheostomy Tubes have a radiopaque, biocompatible outer cannula constructed of polyvinyl chloride. A swivel neck plate allows conformity to individual neck anatomies. Shiley Tracheostomy Tubes (LPC, FEN, CFS, CFN) are available in four sizes: 4, 6, 8, and 10. The Shiley Laryngectomy Tube (LGT) is available in three sizes: 6, 8, and 10. The reusable inner cannula with integral white 15 mm twist-lock connector can be used for standard respiratory equipment and is translucent for easy inspection. The smooth, rounded tip obturator facilitates insertion.

Cuffed models (LPC, FEN) have a thin wall, high volume, low pressure cuff to minimize pressure in the trachea. When inflated, the cuff conforms to the natural shape of the trachea providing a seal at low intracuff pressure. The cuff inflation line has a luer valve with an integral pilot balloon to indicate cuff inflation.

Fenestrated models (CFN, FEN) also feature a white 15 mm Cap (CAP) and a reusable fenestrated inner cannula with green 15 mm twist-lock connector. These accessories can be used to facilitate phonation with or without ancillary speaking devices or for use during the weaning process. The fenestrated inner cannula with green 15 mm connector is not for use when mechanical ventilation is required.

Cuffless models (CFS, CFN, LGT) also feature inner cannulae with low-profile connectors. This type of inner cannula is designed to reduce tube protrusion at the neck but is not compatible with respiratory equipment.

The red Decannulation Plug (DCP) can be used to occlude the proximal end of the outer cannula, forcing the patient to breathe through the fenestrations and upper airway tract during the weaning process. The Decannulation Plug is available in all four sizes of the Shiley Tracheostomy Tubes and can be purchased separately.

#### **Indications for Use**

The devices are intended for use in providing tracheal access for airway management. The fenestrated devices (FEN, CFN) are also indicated when the use of fenestrations is desirable in order to safely and effectively wean a tracheostomy patient. When used in conjunction with the Decannulation Plug (DCP), the fenestrated tracheostomy tube can provide a means of weaning and/or phonation for the patient. Use of the Decannulation Plug (DCP) with the fenestrated tracheostomy tubes forces air through the fenestrations and around the tube, and into the upper airway and vocal cords.

#### **CAUTION:**

*The Shiley™ Tracheostomy Tube is classified as a disposable medical device. The manufacturer recommends that the Tracheostomy Tube usage not exceed twenty-nine (29) days. Frequent and routine changes of the Tracheostomy Tube and accessories are recommended and should be evaluated by the attending physician.*

#### **Instructions**

##### **Tube Preparation**

1. Selection of the tracheostomy tube size is left to the discretion of the physician. Patients in the home care environment should be carefully instructed by their home health care provider in the proper use and handling of this device.

##### **Pre-Insertion Cuff and Inflation Test**

**NOTE:** Refer to Table 1 for Leak Test Inflation Volumes. Inflation volumes are for test purposes only. Consult the physician or Home Health Care Provider for the appropriate inflation volume/pressure when the tube is positioned within the trachea.

2. With Shiley Cuffed models (LPC, FEN) the cuff and inflation system should be tested for leakage before inserting the tube. This test can be performed as follows: Inflate the cuff with the volume of air indicated in Table 1. Then either observe for deflation over several minutes or immerse the tube in sterile saline and observe for air leakage. Deflate the cuff prior to insertion.

## Insertion

### **CAUTION:**

*To ease insertion and to guard against cuff perforation from sharp edges of cartilage, the cuff should be tapered back. This can be accomplished by first inflating the cuff. Then gently move the cuff away from the distal tip of the outer cannula towards the swivel neck plate as the residual air is removed by deflation. Do not use any sharp instruments such as forceps or hemostats that would damage the cuff when tapering it.*

3. Remove the non-fenestrated inner cannula with the white 15mm connector and insert the obturator into the outer cannula. The obturator should be fully seated before the tube is inserted into the patient. A thin film of water soluble lubricant can be applied to the outer cannula, cuff (LPC, FEN), and protruding portion of the obturator to facilitate insertion.
4. After the tracheotomy procedure is performed, insert the tube into the patient's trachea. After the tube is properly in place, remove the obturator immediately.

### **CAUTIONS:**

*When using a water soluble lubricant verify that the lubricant does not enter and occlude the tube lumen which may restrict ventilation.*

*The obturator should be cleaned and then stored in an accessible location near the patient should it be needed for an unscheduled re-intubation.*

### **WARNINGS:**

**ALL INNER CANNULAE WITH TWIST-LOCK CONNECTORS PROVIDED IN THIS PACKAGE SHOULD NOT BE USED ON ANY OTHER TUBE BECAUSE THEY ARE FITTED TO THE EXACT LENGTH OF THIS TUBE.**

**DURING AND AFTER ATTACHMENT OF RESPIRATORY OR ANESTHESIA TUBING AND/OR CONNECTORS TO THE INNER CANNULA, AVOID APPLICATION OF EXCESSIVE ROTATIONAL, LINEAR, OR ROCKING FORCES ON THE TUBING AND/OR CONNECTORS TO PREVENT ACCIDENTAL DISCONNECTION OF THE INNER CANNULA, OR DAMAGE TO THE TRACHEOSTOMY TUBE.**

5. Prior to inserting the inner cannula, establish the patency of the patient's upper airway track. The patient's airway should be cleared by coughing and/or suctioning. Moisten the non-fenestrated inner cannula with white 15 mm connector with sterile saline to facilitate insertion. To lock the non-fenestrated inner cannula in place, hold the swivel neck plate firmly with your fingertips and twist the connector clockwise, one-quarter turn past the locking detent. The outer cannula should always be used with the inner cannula in place unless the Decannulation Plug (DCP) is being used.  
**or**
6. To insert the fenestrated inner cannula with the green 15 mm connector (FEN, CFN), establish the patency of the patient's upper airway track. The patient's airway

should be cleared by coughing and/or suctioning before inserting the fenestrated inner cannula. Moisten the fenestrated inner cannula with the green 15 mm connector with sterile saline to facilitate insertion. To lock the fenestrated inner cannula in place, hold the swivel neck plate firmly with your fingertips and twist the green 15 mm connector clockwise one-quarter turn past the locking detent.

7. The patient's respirations and vital signs should be routinely evaluated after locking the inner cannula in the tube.

**CAUTION:**

*To avoid applying pressure against the patient, the neckplate may be stabilized with your free hand during the locking procedure.*

*Verify that the twist-lock connector engages securely after each use. If parts become worn or loose, immediately report this to your physician for prompt replacement of the tracheostomy tube.*

**WARNINGS:**

**DO NOT ATTEMPT TO VENTILATE THE PATIENT WITH THE GREEN 15 mm CONNECTOR FENESTRATED INNER CANNULA IN PLACE.**

**IF YOU ARE UNABLE TO REMOVE THE INNER CANNULA, CONTACT A PHYSICIAN IMMEDIATELY.**

**Cuff Inflation**

**WARNINGS:**

**UNDER NO CIRCUMSTANCES SHOULD MORE THAN 25 mm OF MERCURY AIR PRESSURE BE USED TO INFLATE THE CUFF. OVER-INFLATION OF THE CUFF MAY CAUSE TRACHEAL DAMAGE AND MAY INHIBIT VENTILATION.**

**AS A FURTHER PRECAUTION FOR VENTILATOR-DEPENDENT PATIENTS, CUFF INFLATION SHOULD BE CHECKED ON A REGULAR BASIS AND REPLACEMENT TRACHEOSTOMY TUBES SHOULD BE KEPT AT BEDSIDE.**

8. Inflate the low pressure cuff by injecting air into the luer valve of the inflation line using a syringe. Selection of a cuff inflation and deflation procedure is left to the discretion of the physician.

**CAUTION:**

*These cuffed products (LPC, FEN) are composed of soft materials to conform to tracheal tissue for performance and patient comfort. Simple precautions in handling of the Shiley™ Cuffed Tracheostomy Tubes during insertion and while in place will facilitate proper function and minimize tears and breaks in the inflation system. Avoid pulling or manipulation of the inflation line, as it is designed to conduct and hold air as part of the cuff inflation system. It is recommended that the inflation line be maintained in a position allowing for patient mobility without placing tension on the line-to-cannula junction. Prevent lint or other particulates from entering the luer valve of the pilot balloon.*

9. Secure the tracheostomy tube to the patient using the neck strap provided.

#### Cuff Deflation

10. Accumulated secretions above the cuff (LPC, FEN) may need to be suctioned before deflating the cuff with a syringe, unless suctioning is contraindicated.
11. To deflate the low pressure cuff (LPC, FEN), withdraw air slowly from the luer valve of the inflation line using a syringe.

#### Cleaning

##### WARNINGS:

**DO NOT USE SOLUTIONS OR CHEMICAL AGENTS OTHER THAN THOSE RECOMMENDED IN THE TABLE BELOW TO CLEAN ANY PART OF THE TRACHEOSTOMY TUBE, AS THIS MAY RESULT IN TUBE DAMAGE.**

**DO NOT SOAK ANY PART OF THE TUBE IN HYDROGEN PEROXIDE OR ANY OTHER SOLUTION.**

• CUFFED OUTER CANNULA	RINSE GENTLY WITH STERILE SALINE
• CUFFLESS OUTER CANNULA	1. CLEAN WITH HYDROGEN PEROXIDE (HALF-STRENGTH), STERILE NORMAL SALINE, OR WATER AND MILD DETERGENT.
• NECK PLATE	2. AFTER CLEANING THOROUGHLY RINSE WITH STERILE SALINE TO REMOVE ANY CLEANING SOLUTION RESIDUES.
• INNER CANNULA	3. AIR DRY.
• OBTURATOR	
• CAP	
• DECANNULATION PLUG	

#### CAP (FEN, CFN)

The white 15 mm Cap is for use with the fenestrated inner cannula with the green 15 mm connector. The CAP when attached to the proximal end of the green 15 mm connector of the fenestrated inner cannula will block all expired and inspired air through the tube, forcing breathing through the mouth and nose, allowing speech.

#### Decannulation Plug (DCP)

The Decannulation Plug (DCP) is red for easy identification and is available in all four sizes 4, 6, 8 and 10 of the Shiley™ Tracheostomy Tube models. The Decannulation Plug occludes the proximal end of the outer cannula, forcing the patient to breathe through the fenestrations and upper airway tract during the weaning process.

**WARNINGS:**

**USE THE DECANNULATION PLUG (DCP) AND CAP ONLY WITH SHILEY FENESTRATED LOW PRESSURE CUFFED (FEN), AND CUFFLESS FENESTRATED (CFN) TRACHEOSTOMY TUBES.**

**PRIOR TO INSERTING THE DECANNULATION PLUG (DCP) OR ATTACHING THE 15 mm CAP TO THE FENESTRATED INNER CANNULA WITH THE GREEN 15 mm CONNECTOR, MAKE SURE THAT THE FENESTRATIONS ARE NOT OCCLUDED, THAT THE CUFF (FEN) IS COMPLETELY DEFLATED AND THAT THERE IS SUFFICIENT AIRWAY FOR THE PATIENT.**

12. Establish the patency of the patient's upper airway tract. The patient's airway should be cleared by coughing and/or suctioning before plugging the tracheostomy tube with the Decannulation Plug or capping the fenestrated inner cannula with the green 15 mm connector to the 15 mm Cap (CAP).
13. Completely deflate the tracheostomy tube cuff (FEN only).
14. To use the 15 mm Cap (CAP):  
Attach the 15 mm Cap to the green 15mm connector of the fenestrated inner cannula.

**or**

15. To use the DCP:  
Remove the inner cannula with 15mm twist-lock connector and insert the Decannulation Plug with twist-lock connector to the proximal end of the outer cannula.
16. To lock the Decannulation Plug in place, hold the swivel neck plate firmly with your fingertips and twist the Decannulation Plug clockwise past the locking detent one quarter turn.

***CAUTION: Locking the Decannulation Plug on an angle could strip the locking mechanism.***

17. The patient's respirations and vital signs should be evaluated after plugging or capping the tracheostomy tube.
18. If the patient needs to be mechanically ventilated while using the Decannulation Plug or the 15 mm Cap, replace the accessory with the non-fenestrated inner cannula with the white 15 mm twist-lock connector, inflate the cuff (FEN), and connect the tracheostomy tube to the ventilator.

**WARNING:**

**IF THE PATIENT HAS DIFFICULTY BREATHING WHILE USING THE DECANNULATION PLUG OR THE 15 mm CAP / FENESTRATED INNER CANNULA WITH GREEN 15 mm CONNECTOR, REMOVE THOSE ACCESSORIES, REPLACE WITH THE NON-FENESTRATED INNER CANNULA IF NECESSARY AND VERIFY AIRWAY PATENCY.**

## Removing The Tracheostomy Tube

19. The inner cannula can be cleaned and reused to maintain a patent airway without requiring the entire tube to be changed.

### **WARNINGS:**

**SPARE INNER CANNULAE (SIC) WITH RED 15 mm CONNECTORS ARE AVAILABLE SEPARATELY FOR SHORT TERM USE TO MAINTAIN MECHANICAL VENTILATION WHILE THE ORIGINAL INNER CANNULA IS BEING CLEANED. IT SHOULD NOT BE USED AS A REPLACEMENT FOR THE INNER CANNULA PROVIDED WITH THE PRODUCT.**

20. Before removing the tracheostomy tube, the cuff (if present) should be completely deflated with a syringe. This will ensure that the cuff passes through the stoma with minimum resistance.

## Additional Instructions for Use - Fenestrated Tracheostomy Tubes (FEN, CFN)

### **WARNINGS:**

**SPECIAL CARE MUST BE EXERCISED TO MAINTAIN FENESTRATION PATENCY, PARTICULARLY FOR PATIENTS WITH THICK, TENACIOUS SECRETIONS THAT MAY BLOCK THE FENESTRATIONS. USE APPROPRIATE HUMIDIFICATION WHEN INDICATED TO ASSIST IN MAINTAINING LUMEN PATENCY. CARE SHOULD BE TAKEN IN THE SELECTION AND PLACEMENT OF THE FEN OR CFN TRACHEOSTOMY TUBE TO MINIMIZE THE COMPLICATIONS THAT MAY ARISE FROM GRANULATION TISSUE. OBSTRUCTION OF THE FENESTRATIONS MAY PREVENT INSERTION OF THE NON-FENESTRATED INNER CANNULA WITH THE WHITE 15 mm CONNECTOR NEEDED FOR MECHANICAL VENTILATION OR COMPROMISE UPPER AIRWAY ACCESS. IF THIS OCCURS, THE FEN OR CFN TRACHEOSTOMY TUBE MAY NEED TO BE REPLACED WITH ANOTHER TYPE OF TUBE.**

**WHEN USING A FENESTRATED TUBE (FEN or CFN) THE NON-FENESTRATED INNER CANNULA WITH THE WHITE 15 mm CONNECTOR MUST BE USED DURING SUCTIONING TO PREVENT THE SUCTION CATHETER FROM PROTRUDING THROUGH THE FENESTRATIONS OF THE TRACHEOSTOMY TUBE. SUCTIONING WITHOUT THE NON-FENESTRATED INNER CANNULA IN PLACE COULD RESULT IN DAMAGE TO THE TRACHEAL WALL AND COULD CAUSE THE SUCTION CATHETER TO BECOME CAUGHT IN THE FENESTRATIONS. IF YOU FEEL RESISTANCE WHEN INSERTING THE INNER CANNULA PAST THE FENESTRATION, DO NOT FORCE THE INNER CANNULA THROUGH THE CANNULA. NOTIFY YOUR HEALTH CARE PROVIDER IMMEDIATELY.**

**NEVER USE THE FENESTRATED TRACHEOSTOMY TUBES FOR MECHANICAL VENTILATION UNLESS THE NON FENESTRATED INNER CANNULA WITH WHITE 15 mm CONNECTOR IS IN PLACE TO PREVENT LEAKAGE THROUGH THE FENESTRATIONS.**

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© 2011 Covidien.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
**EC REP** Covidien Ireland Limited,