



# Shiley™

## Adult Flexible Tracheostomy Tube with TaperGuard™ Cuff Disposable Inner Cannula

|           |                          |                    |
|-----------|--------------------------|--------------------|
| <b>en</b> | Instructions For Use     | <b>REF</b> 4CN65A  |
| <b>fr</b> | Mode d'emploi            | <b>REF</b> 5CN70A  |
| <b>de</b> | Gebrauchsanleitung       | <b>REF</b> 6CN75A  |
| <b>nl</b> | Gebruiksaanwijzing       | <b>REF</b> 7CN80A  |
| <b>it</b> | Istruzioni per l'uso     | <b>REF</b> 8CN85A  |
| <b>es</b> | Instrucciones de uso     | <b>REF</b> 9CN90A  |
| <b>sv</b> | Bruksanvisning           | <b>REF</b> 10CN10A |
| <b>da</b> | Brugsanvisning           |                    |
| <b>no</b> | Bruksanvisning           |                    |
| <b>fi</b> | Käyttöohjeet             |                    |
| <b>pt</b> | Instruções de uso        |                    |
| <b>ru</b> | Инструкции по применению |                    |
| <b>pl</b> | Instrukcja użytkowania   |                    |
| <b>sk</b> | Návod na použitie        |                    |
| <b>sl</b> | Navodilo za uporabo      |                    |
| <b>el</b> | Οδηγίες χρήσης           |                    |
| <b>tr</b> | Kullanım Talimatları     |                    |
| <b>lt</b> | Naudojimo instrukcija    |                    |

| REF     |      |      |    | C (mm) |    | (cc) |  | (mm) |    | (mm) |  |  |  |  |  |
|---------|------|------|----|--------|----|------|--|------|----|------|--|--|--|--|--|
|         |      |      |    |        |    |      |  |      |    |      |  |  |  |  |  |
| 4CN65A  | 6.5  | 9.4  | 73 | 62     | 11 | 20.6 |  | 5.5  | 62 |      |  |  |  |  |  |
| 5CN70A  | 7.0  | 10.1 | 79 | 68     | 13 | 23.0 |  | 6.0  | 68 |      |  |  |  |  |  |
| 6CN75A  | 7.5  | 10.8 | 85 | 74     | 14 | 25.4 |  | 6.5  | 74 |      |  |  |  |  |  |
| 7CN80A  | 8.0  | 11.4 | 88 | 77     | 16 | 25.4 |  | 7.0  | 77 |      |  |  |  |  |  |
| 8CN85A  | 8.5  | 12.2 | 90 | 79     | 17 | 26.6 |  | 7.5  | 79 |      |  |  |  |  |  |
| 9CN90A  | 9.0  | 12.7 | 90 | 79     | 19 | 27.6 |  | 8.0  | 79 |      |  |  |  |  |  |
| 10CNT0A | 10.0 | 13.8 | 90 | 79     | 20 | 28.6 |  | 9.0  | 79 |      |  |  |  |  |  |

| REF |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  | Part Number, Numéro de la pièce, Einzelteilnummer, Onderdeelnummer, Codice, Número de peça, Artikelnummer, Delnummer, Delenummer, Osanummer, Número da peça, Номер по каталогу, Numer katalogowy, Císlo dielü, Številka dela, Αριθμός προϊόντος, Παρçá Numarasi, Dalles   |
|-----|--|--|--|--|--|--|--|--|--|--|--|--|--|--|--|---|
|     |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  | Part Number, Numéro de la pièce, Einzelteilnummer, Onderdeelnummer, Codice, Número de peça, Artikelnummer, Delnummer, Delenummer, Osanummer, Número da peça, Номер по каталогу, Numer katalogowy, Císlo dielü, Številka dela, Αριθμός προϊόντος, Παρçá Numarasi, Dalles   |
|     |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  | ID – Inside Diameter (mm), DI – Diamètre interne (mm), ID – Innendurchmesser (mm), ID – Innendurchmesser (mm), ID – binnendiameter (mm), ID – Diametro interno (mm), OD – Diametro estero (mm), DE – Diametro exterior (mm), YD – Ytterdiameter (mm), ID – indvending diameter (mm), ID – indvending diameter (mm), ID – Diámetro interno (mm), DI – Diámetro interno (mm), ВД – внутренний диаметр (mm), ID – średnica wewnętrzna (mm), ID – średnica wewnętrzna (mm), ID – Vnitorný priemer (mm), ID – notranji premer (mm), ID – Eσωτερική διάμετρος (mm), ID – İç Çap (mm), ID – vidinis diametras (mm) |
|     |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  | OD – Outside Diameter (mm), DE – Diamètre externe (mm), AD – Außendurchmesser (mm), OD – buitendiameter (mm), OD – Diametro estero (mm), DE – Diametro exterior (mm), YD – Ytterdiameter (mm), OD – udvending diameter (mm), OD – Diámetro externo (mm), OD – ulkohalkaisija (mm), DE – Diámetro Externo (mm), HD – наружный диаметр (mm), OD – średnica zewnętrzna (mm), OD – Vonkajši premer (mm), OD – zunanji premer (mm), OD – εξωτερική διάμετρος (mm), OD – Diş Çap (mm), OD – išorinis diametras (mm)   |

|   |  |
|---|--|
|   | <p>C – Tube Length (mm), C – Longueur du tube (mm), C – Kanülenlänge (mm), C – tubelengte (mm), C – Lunghezza cannula (mm), C – Longitud del tubo (mm), C – Tublängd (mm), C – tubens længde (mm), C – Tubelengde (mm), C – putken pituus (mm), C – Comprimento do tubo (mm), C – длина трубки (мм), C – dlugosc rurki (mm), C – Džška trubice (mm), C – dolžina cevi (mm), C – Μηκος σωλήνα (mm), C – Túp Uzunluđu (mm), C – vamzdeljo ilgis (mm)</p> <p><math>\theta</math> – Tube Angle (°), <math>\theta</math> – Angle du tube (°), <math>\theta</math> – Krümmung der Kanüle (°), <math>\theta</math> – tube-hoek (°), <math>\theta</math> – Angulo del tubo (°), <math>\theta</math> – Tubens vinkel (°), <math>\theta</math> – tubens vinkel (°), <math>\theta</math> – putken kulma (°), <math>\theta</math> – Angulo do tubo (°), <math>\theta</math> – угол трубки (°), <math>\theta</math> – kät rurki (°), <math>\theta</math> – Zakrivnenje trubice (°), <math>\theta</math> – kot cevi (°), <math>\theta</math> – Γούια σωλήνα (°), <math>\theta</math> – Túp Açısı (°), <math>\theta</math> – vamzdelio kampas (°)</p> |
|  | <p>Disposable Inner Cannula Length (mm), Longueur de la canule interne jetable (mm), Länge Innenkanüle zum Einmalgebrauch (mm), Länge wegwerbare Binnencanüle (mm), Lunghezza contro cannula monouso (mm), Longitud de la cánula interna desechable (mm), Längd på innekanyli för engångsbruk (mm), Længden af inderkanylen til engangsbrug (mm), Engangsinnerekanyle lengde (mm), Kertäkäyttöisen sisäkanyylin pituus (mm), Comprimento da canula interna descartável (mm), Длина одноразовой внутренней канюли (мм), Długość wewnętrznej kaniumi jednorazowego użytku (mm), Džška jednorazovej kanule za enkratno uporabo (mm), Μηκος εσωτερικής κανούλας μιας χρήσης (mm), Tek kullanimilik iç Kanül Uzunluđu (mm), Vienkartinis vidlines kaniume ilgis (mm)</p>  |
|  | <p>Test Inflation Volume (cc), syringe not included, Volume de gonflage de test (cc) ; seringue non fournie, Test-Befüllungsvolumen (cc), Spritze nicht enthalten, Inflationevolume test (cc); spuit niet inbegrepen, Volume di gonfiaggio per il test (cc); siringa non inclusa, Volumen de inflado de prueba (cc); jeringa no incluida, Volym för testupblåsning (ml), spruta ingår inte, Testinflationsvolumen (cc), sprøjte medfølger ikke, Test fyllvolum (cc), sprøjte ikke inkludert, Testitältötilavuus (cc), ei sisällytä ruiskua, Volume do teste de insuflação (cc); seringa não incluida, Testovaci objem nakuchivaniya (см3); шприц не входит в комплект, Objętość napełnienia testowego (ml); brak strzykawki w zestawie, Testovací objem nafúknutia (cc); striekačka nie je súčasťou balenia, Prostornina za preskusno naphovanje (cc); injekcija ni vključena, Докиматікс оукоқ дақстолық (cc); Δεν περιλαμβάνεται σύριγγα, Şişirme Testi Hacmi (cc); şiริงga dahil değildir, Bandomo pūrimo tūris (cc); svirkštas nepridedamas</p>  |
|  | <p>Cuff Resting Diameter (mm), Diamètre du ballonnet au repos (mm), Cuff-Durchmesser, entleert (mm), Nominale diameter cuff (mm), Diametro residuo della cuffia (mm), Diámetro nominal del balón (mm), Kuffens diameter i vila (mm), Diameter på deflateret cuff (mm), Mansjett hvilende diameter (mm), Kuffin lepoalkaisija (mm), Diámetro do apolo do balão (mm), Диаметр пустой манжеты (мм), Średnica spoczynkowa mankietu (mm), Kludový priemer manžety (mm), Premier nenapuhnjene manšete (mm), Διαμέτρος αεροβάλλου σε κατάσταση ηρεμίας (mm), Kaf Oturma Çapı (mm), Nepripūstos manžetės diametras (mm)</p>  |
|  | <p>Loading Dilator (fr); not included, Dilatateur de charge (fr); non fourni, Ladedilatator (fr); nicht inbegriffen, Laaddilatator (fr), niet inbegrepen, Introductore di cateteri a palloncino (fr); non incluso, Dilatador de carga (fr); no incluido, Dilator (fr), ingår inte, Indföringsdilatator (fr), medfølger ikke, Fylldilatator (fr), ikke inkludert, Asetusdilatattori (fr), ei kuulu mukaan, Dilatador de carregamento (fr); não incluido, Дилататор (fr); не входит в комплект, Rozszerzadło (fr), brak w zestawie, Dilatátor zavádzania (fr); nie je súčasťou balenia, Pohlinli dilatator (fr); ni priložen, Διαστολέας τοποθέτησης (fr); Δεν περιλαμβάνεται, Yükleme Dilatörü (fr); dahil değildir, Jvediklis (fr); nepridedamas</p>   |
|  | <p>Obturator, Introducteur, Obturator, Obturatore, Obturador, Obturator, Obturador, Obturador, Obturaattori, Obturator, Obturaattori, Obturador, Обтуратор, Obturator, Obturateur, Testnilo, Επιμοιωτικό, Obturator, Kaišis</p>  |
|  | <p>Neck Strap, Attache cervicale, Halteband, Nekkband, Fettuccia di fissaggio, Cinta para el cuello, Nackband, Halsrem, Halsstropp, Niskanauha, Correia para pescoço, Лента для крепления на шею, Pasek szyjny, Popruha na krk, Trak za pritrditev na vrat, Λωρίδα υποστηρίξης αυχένα, Bounyn Kayış, Kaklo raštis</p>  |

# Shiley™

## Adult Flexible Tracheostomy Tube with TaperGuard™ Cuff

### Disposable Inner Cannula

Identification of a substance that is contained or present within the product or packaging.

Identification of a substance that is not contained or present within the product or packaging.

|            |         |         |
|------------|---------|---------|
| <b>REF</b> | 4CN65A  | 6.5 mm  |
| <b>REF</b> | 5CN70A  | 7.0 mm  |
| <b>REF</b> | 6CN75A  | 7.5 mm  |
| <b>REF</b> | 7CN80A  | 8.0 mm  |
| <b>REF</b> | 8CN85A  | 8.5 mm  |
| <b>REF</b> | 9CN90A  | 9.0 mm  |
| <b>REF</b> | 10CN10A | 10.0 mm |

### Instructions for Use

POST FOR USE BY ALL TRACHEOSTOMY CARE PERSONNEL

#### Description

The Shiley™ Tracheostomy Tube with TaperGuard™ Cuff (“tracheostomy tube”) is a dual cannula tracheostomy tube with a disposable inner cannula. It contains an integrated 15 mm connector for use with standard ventilation and anesthesia equipment.

The tracheostomy tube has a radiopaque, biocompatible outer cannula constructed of polyvinyl chloride. A pliable neck plate allows conformity to individual neck anatomies and contains two suture holes. The tapered tracheostomy tube tip is compatible with percutaneous insertion. The smooth, rounded tip obturator facilitates insertion. The disposable inner cannula is translucent for easy inspection and allows for airway maintenance.

The tracheostomy tube has a taper-shaped, low-pressure cuff to provide an effective air and fluid seal and minimize tracheal wall damage. The cuff inflation line has a luer valve with an integral pilot balloon.

This tracheostomy tube can be used with or without the disposable inner cannula included in the package (and also available separately; reference 10006780 for part numbers). However, use of the inner cannula is recommended.

#### Indications for Use

The Shiley™ Tracheostomy Tube with TaperGuard™ Cuff is intended to provide tracheal access for airway management.

The Shiley™ Tracheostomy Tube with TaperGuard™ Cuff is also intended for use with Percutaneous Dilatational Tracheotomy (PDT) procedures.

#### WARNINGS

- Under no circumstances should more than 25 cm of H<sub>2</sub>O air pressure be used to inflate the cuff. Over-inflation of the cuff may cause tracheal damage and may inhibit ventilation. Consult with the attending physician.
- Test each tube’s cuff, pilot balloon, and valve by inflation prior to use. If dysfunction is detected in any part of the inflation system, the tube should not be used.
- Inflation volumes listed on page 2 are for test purposes only. Consult the physician or home health care provider for the appropriate inflation volume/pressure when the tube is positioned within the trachea.
- When using a laser or electrosurgical device in conjunction with this or other tracheal tubes, care must be taken to avoid contact with the tube. Such contact, especially in the presence of oxygen-enriched mixtures or nitrous oxide containing mixtures, could result in rapid combustion of the tube with harmful thermal effects and with emission of corrosive and toxic combustion products, including hydrochloric acid (HCl).
- Diffusion of nitrous oxide, oxygen, or air can alter cuff volume and pressure.
- The tracheostomy tube and neck strap cannot be adequately resterilized by the user in order to facilitate safe reuse and are therefore intended for single patient use. Attempts to resterilize

this device may result in product failure and increased risks to the patient.

- Sterile only if the protective tray and cover are not open, damaged, or broken. Do not resterilize. Do not expose to temperatures above 120°F (49°C).
- Prior to use, handle with due care to prevent user contamination.
- The disposable inner cannula is designed for single use and should not be cleaned or reused.
- When using a disposable inner cannula, it can only be replaced by a Shiley disposable inner cannula of the same size.
- It is recommended that a spare disposable inner cannula be with the patient at all times.
- When using a water-soluble lubricant, verify that the lubricant does not enter and occlude the tube lumen, which may restrict ventilation.
- As a further precaution for ventilator-dependent patients, cuff inflation should be checked on a regular basis, and replacement tracheostomy tubes should be kept at bedside.
- Syringes, three-way stopcocks, or other devices should not be left inserted in the inflation valve for an extended period of time. The resulting stress could crack the valve housing and allow the cuff to deflate.

### CAUTIONS

- Federal law (USA) restricts this device to sale by or on the order of a physician.
- The tracheostomy tube and obturator are single patient use medical devices. Duration of single patient use should not exceed twenty-nine (29) days. Covidien does not recommend and has not substantiated use of these devices beyond the 29-day time frame. Decisions about tracheostomy tube changes should be made by the responsible physician or designate using accepted medical techniques and judgment.
- Clinicians placing this device percutaneously should be competent in the procedure prior to using this tube.
- Prior to use, check that the obturator can be inserted and removed without difficulty.
- To ease placement of the tracheostomy tube, generously lubricate the surface of the loading dilator.
- To ease insertion and to guard against cuff perforation from sharp edges of cartilage, taper back the cuff by deflating the cuff and gently moving it away from the distal tip of the cannula toward the neck plate as the residual air is removed by deflation. When tapering the cuff, do not use sharp instruments that would damage the cuff, such as forceps or hemostats.
- Avoid pulling or manipulation of the inflation line, as it is designed to conduct and hold air as part of the cuff inflation system.
- Prevent lint or other particulates from entering the luer valve of the pilot balloon.
- Extremes of movement may result in mucosal abrasion and transient disruption of cuff seal.
- The use of Lidocaine Topical Aerosol has been associated with the formation of pin holes in PVC cuffs. Lidocaine hydrochloride solution has been reported not to have this effect.
- Do not use solutions or chemical agents other than those recommended in these instructions to clean any part of the tracheostomy tube, as this may result in tube damage.
- Precautions should be taken when discarding the tracheostomy tube and accessories. Disposal should be made in accordance with current institutional standards and applicable national regulations for biologically hazardous waste.

### NOTES

- The tracheostomy tube and accessories are safe to use for X-ray, ultrasound, Positron Emission Tomography (PET), and Radiation Therapy.
- During an MRI scan, the pilot balloon should be secured near the Y connector of the ventilator circuit at least 3 cm from the area of interest to prevent movement and image distortion.
- The tracheostomy tube and accessories are composed of medical grade plastics which do not contain and have not been manufactured with natural rubber latex material.
- Sterilized with ethylene oxide.

**Contents of Package:** 1 tracheostomy tube with cuff and pliable neck plate; 1 disposable inner can-

nula; 1 obturator; 1 neck strap. Reference page 2 for sizes and dimensions.

**Accessories Available from Covidien:** CAP, Disposable Inner Cannula (reference 10006780 for part numbers).

**Accessories Available Separately:** Some of the tube sizes (reference page 2) are compatible with the Percutaneous Dilatational Loading Dilator from the COOK® Critical Care Ciaglia Percutaneous Introducer Sets and Trays.

**NOTE:** The sizing dimensions listed on page 2 refer to the ID (the inside diameter of the tube or inner cannula at its narrowest point), the OD (the outside diameter of the tube, not including cuff features), and the length (the distance from the neck plate to the distal tip on the tube's or inner cannula's centerline).

### Tube Preparation and Cuff Inflation Test

1. Based on the physician's discretion, select the appropriate tracheostomy tube size.
2. Test the cuff and inflation system for leakage before inserting the tube, as follows:
  - a. Inflate the cuff with the volume of air indicated on page 2.
  - b. Either observe for deflation over several minutes, or immerse the tube in sterile saline and observe for air leakage.
  - c. Deflate the cuff before insertion.

### Tube Insertion Using the Obturator or a Percutaneous Loading Dilator

The tube sizes that have a compatible loading dilator available (reference page 2) can be used in percutaneous dilatational tracheotomy (PDT) procedures.

**NOTE:** As this tracheostomy tube can be placed percutaneously, sterile technique should be followed for the handling and placement of the tracheostomy tube.

1. To ease insertion and to guard against cuff perforation from sharp edges of cartilage, taper back the cuff as follows:
  - a. Deflate the cuff.
  - b. Gently move the cuff away from the distal tip of the tracheostomy tube toward the neck plate as the residual air is removed by deflation.
2. Insert the **obturator or compatible loading dilator** into the tracheostomy tube.
  - If using the **obturator**, ensure it is fully seated. If indicated to ease insertion, apply a thin film of water soluble lubricant to the outside of the tracheostomy tube, cuff, and protruding portion of the obturator.
  - If using a **loading dilator for a PDT procedure**, generously lubricate the surface of the appropriately sized loading dilator and load the tracheostomy tube onto the dilator. Ensure that the tapered section of the loading dilator clears the distal tip of the tracheostomy tube by approximately 2 cm.

**NOTE:** The loading dilator is available separately.

3. If using a guide wire with the obturator, first insert the guide wire in the patient, then thread it all the way through the obturator, starting at the distal tip of the obturator.

### WARNING

If the guide wire is damaged during the procedure, it may be difficult to continue insertion without trauma to the trachea and misplacement of the tracheostomy tube. In these cases, it may be possible to continue only if:

- a. the damaged section can be advanced into the trachea, and
- b. there is a sufficient length of undamaged guide wire remaining.

If the guide wire cannot be salvaged, use a new guide wire.

4. Insert the tube into the patient's trachea.
5. After the tube is properly in place, remove the **loading dilator or obturator** and guide wire (if used) immediately.
6. Insert the disposable inner cannula into the tracheostomy tube and press to lock into position.
7. Clean the obturator and store it in an accessible location near the patient in case it is needed for an unscheduled re-insertion.

### Cuff Inflation

1. Inflate the cuff only with enough gas mixture to provide an effective seal at the desired lung inflation pressure. Use Minimal Occluding Volume or Minimum Leak techniques in conjunction with an intracuff pressure measuring device in selecting the sealing pressure.
2. Secure the tracheostomy tube to the patient using the neck strap provided.

## Cuff Deflation

To deflate the cuff, use a syringe to slowly withdraw air from the luer valve of the inflation line.

## Tube Removal

1. Remove the neck strap from around the patient's neck.
2. Withdraw air slowly from the luer valve of the inflation line using a syringe. Completely deflate the cuff. Deflating the cuff ensures that the cuff passes through the stoma with minimum resistance.
3. While holding on to the tracheostomy tube's 15 mm connector, gently twist and pull the breathing circuit connector to disconnect the circuit from the tracheostomy tube.
4. Slowly remove the tracheostomy tube from the patient.

## Cleaning

|   |   |
|---|---|
| <ul style="list-style-type: none"><li>• Disposable inner cannula</li><li>• Neck strap</li></ul> | Do not clean or reuse.  |
| <ul style="list-style-type: none"><li>• Cuffed tracheostomy tube</li></ul>                      | <ol style="list-style-type: none"><li>1. Rinse gently with sterile, normal saline. Do not soak any part of the tube in hydrogen peroxide or any other solution.</li><li>2. Air dry.</li></ol>   |
| <ul style="list-style-type: none"><li>• Neck plate</li><li>• Obturator</li></ul>                | <ol style="list-style-type: none"><li>1. Clean with hydrogen peroxide (half-strength); sterile, normal saline; or water and mild detergent. Do not soak any part of the tube in hydrogen peroxide or any other solution.</li><li>2. After cleaning, thoroughly rinse with sterile saline.</li></ol> |

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© 2013 Covidien.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.
Other brands are trademarks of a Covidien company.
Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.
Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.