# Shiley™

## Adult Flexible Evac Tracheostomy Tube with TaperGuard™ Cuff Disposable Inner Cannula

<b>en</b>	Instructions For Use	<b>REF</b>	4CN65ED
<b>fr</b>	Mode d'emploi	<b>REF</b>	5CN70ED
<b>de</b>	Gebrauchsanweisung	<b>REF</b>	6CN75ED
<b>nl</b>	Gebruiksaanwijzing	<b>REF</b>	7CN80ED
<b>it</b>	Istruzioni per l'uso	<b>REF</b>	8CN85ED
<b>es</b>	Instrucciones de uso	<b>REF</b>	9CN90ED
<b>sv</b>	Brugsanvisning	<b>REF</b>	10CN10ED
<b>da</b>	Brugsanvisning		
<b>no</b>	Brugsanvisning		
<b>fi</b>	Käyttöohjeet		
<b>pt</b>	Instruções de uso		
<b>ru</b>	Инструкция по применению		
<b>zh</b>	使用说明		
<b>pl</b>	Instrukcja użycia		
<b>cs</b>	Návod k použití		
<b>sk</b>	Návod na použitie		
<b>hu</b>	Használati utasítás		
<b>el</b>	Οδηγίες χρήσης		
<b>tr</b>	Kullanım Talimatı		
<b>th</b>	คำแนะนำการใช้งาน		
<b>ko</b>	사용 지침		
<b>id</b>	Petunjuk Penggunaan		
<b>ar</b>	تعليمات الاستخدام		


REF	REF		REF		REF		REF		REF		REF		REF		REF		REF		REF	
	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)	(mm)
4CN65ED	6.5	9.4	10.0	73	62	35	42	11	20.6	5.5	62	21	✓							
5CN70ED	7.0	10.1	10.7	79	68	35	48	13	23.0	6.0	68	21	✓							
6CN75ED	7.5	10.8	11.3	85	74	35	53	14	25.4	6.5	74	24	✓							
7CN80ED	8.0	11.4	11.8	88	77	34	54	16	25.4	7.0	77	26	✓							
8CN85ED	8.5	12.2	12.5	90	79	34	56	17	26.6	7.5	79	27	✓							
9CN90ED	9.0	12.7	13.1	90	79	34	56	19	27.6	8.0	79	n/a	✓							
10CN10ED	10.0	13.8	13.8	90	79	33	55	20	28.6	9.0	79	n/a	✓							

REF	REF	REF	REF	REF	REF	REF	REF	REF	REF	REF	REF	REF	REF	REF	REF	REF	REF	REF	REF	REF

**en:** Part Number **fr:** Référence **de:** Artikelnummer **nl:** Onderdeelnummer **it:** Numero codice **es:** Número de pieza **sv:** Artikelnummer  
**da:** Delnummer **no:** Delnummer **fi:** Osanumero **pt:** Número de peça **ru:** Каталожный номер **zh:** 产品编号 **pl:** Numer katalogowy **cs:** Číslo dílu  
**sk:** Číslo dielu **hu:** Cikkszám **el:** Αριθμός προϊόντος **tr:** Parça Numarası **th:** หมายเลขอะไหล่ **ko:** 부품 번호 **id:** Nomor Komponen **ar:** رقم القطعة

**en:** ID – Inside Diameter (mm) **fr:** DI – Diamètre interne (mm) **de:** ID – Innendurchmesser (mm) **nl:** ID – binnendiameter (mm) **it:** ID – Diametro interno (mm) **es:** DI: diámetro interior (mm) **sv:** ID – Innerdiameter (mm) **da:** ID – Indvendig diameter (mm) **no:** ID – Indre diameter (mm) **fi:** ID – sisäläpimitta (mm) **pt:** ID – Diâmetro interno (mm) **ru:** ID – внутренний диаметр (мм) **zh:** ID – 内径 (mm) **pl:** ID – średnica wewnętrzna (mm) **cs:** ID – Vnitřní průměr (mm) **sk:** ID – Vnútorý priemer (mm) **hu:** ID – belső átmérő (mm) **el:** ID – Εσωτερική διάμετρος (mm) **tr:** ID – İç Çap (mm) **th:** ID – เส้นผ่าศูนย์กลางภายใน (มม.) **ko:** ID – 내경 (mm) **id:** ID – Diameter Dalam (mm) **ar:** ID – قطر الداخلي (ملم)

	<p><b>en:</b> OD – Outside Diameter (mm) <b>fr:</b> DE – Diamètre externe (mm) <b>de:</b> OD – Außendurchmesser (mm) <b>nl:</b> OD – buitendiameter (mm) <b>it:</b> OD – Diametro esterno (mm) <b>es:</b> DE: diámetro exterior (mm) <b>sv:</b> YD – Ytterdiameter (mm) <b>da:</b> OD – Udvendig diameter (mm) <b>no:</b> YD – Ytre diameter (mm) <b>fi:</b> OD – ulkoläpmitä (mm) <b>pt:</b> DE – Diâmetro externo (mm) <b>ru:</b> OD — наружный диаметр (мм) <b>zh:</b> OD – 外径 (mm) <b>pl:</b> OD – średnica zewnętrzna (mm) <b>cs:</b> OD – Vnější průměr (mm) <b>sk:</b> OD – Vonkajší priemer (mm) <b>hu:</b> OD – külső átmérő (mm) <b>el:</b> OD – Εξωτερική διάμετρος (mm) <b>tr:</b> OD – Dış Çap (mm) <b>th:</b> OD – เส้นผ่าศูนย์กลาง (ม.) <b>ko:</b> OD – 외경 (mm) <b>id:</b> OD – Diameter Luar (mm) <b>ar:</b> OD – قطر الخارجي (مم)</p>
	<p><b>en:</b> OD – Maximum Outside Diameter (mm) <b>fr:</b> DE – Diamètre externe maximal (mm) <b>de:</b> OD – Maximaler Außendurchmesser (mm) <b>nl:</b> OD – maximale buitendiameter (mm) <b>it:</b> OD – Diametro esterno massimo (mm) <b>es:</b> DEM: diámetro exterior máximo (mm) <b>sv:</b> YD – Maximal ytterdiameter (mm) <b>da:</b> OD – Maks. udvendig diameter (mm) <b>no:</b> YD – Maksimum ytre diameter (mm) <b>fi:</b> OD – maksimiulkoläpmitä (mm) <b>pt:</b> DE – Diâmetro externo máximo (mm) <b>ru:</b> OD — максимальный наружный диаметр (мм) <b>zh:</b> OD – 最大外径 (mm) <b>pl:</b> OD – maksymalna średnica zewnętrzna (mm) <b>cs:</b> OD – Maximální vnější průměr (mm) <b>sk:</b> OD – Maximálny vonkajší priemer (mm) <b>hu:</b> OD – maximális külső átmérő (mm) <b>el:</b> OD – Μέγιστη εξωτερική διάμετρος (mm) <b>tr:</b> OD – Maksimum Dış Çap (mm) <b>th:</b> OD – เส้นผ่าศูนย์กลางมากที่สุด (ม.) <b>ko:</b> OD – 최대 외경 (mm) <b>id:</b> OD – Diameter Luar Maksimum (mm) <b>ar:</b> OD – قطر الخارجي الأقصى (مم)</p>
 <p>A = 0 B = 0</p>	<p><b>en:</b> C – Tube Length (mm) <b>fr:</b> C – Longueur du tube (mm) <b>de:</b> C – Kanülenlänge (mm) <b>nl:</b> C – tubelengte (mm) <b>it:</b> C – Lunghezza cannula (mm) <b>es:</b> C: longitud del tubo (mm) <b>sv:</b> C – Tublängd (mm) <b>da:</b> C – Tubens længde (mm) <b>no:</b> C – Tubelengde (mm) <b>fi:</b> C – putken pituus (mm) <b>pt:</b> C – Comprimento do tubo (mm) <b>ru:</b> C — длина трубки (мм) <b>zh:</b> C – 插管长度 (mm) <b>pl:</b> C – długość rurki (mm) <b>cs:</b> C – Délka trubice (mm) <b>sk:</b> C – Dĺžka kanyly (mm) <b>hu:</b> C – tubus hossza (mm) <b>el:</b> C – Μήκος σωλήνα (mm) <b>tr:</b> C – Tüp Uzunluğu (mm) <b>th:</b> C – ทüp Uzunluđu (mm) <b>ko:</b> C – 튜브 길이(mm) <b>id:</b> C – Panjang Tabung (mm) <b>ar:</b> C – طول الأنبوب (مم)</p> <p><b>en:</b> θ – Tube Angle (°) <b>fr:</b> θ – Angle du tube (°) <b>de:</b> θ – Krümmung der Kanüle (°) <b>nl:</b> θ – tubehoek (°) <b>it:</b> θ – Angolo cannula (°) <b>es:</b> θ: ángulo del tubo (°) <b>sv:</b> θ – Tubens vinkel (°) <b>da:</b> θ – Tubens vinkel (°) <b>no:</b> θ – XxxX XxxX (°) <b>fi:</b> θ – putken kulma (°) <b>pt:</b> θ – Angulo do tubo (°) <b>ru:</b> θ — угол изгиба трубки (°) <b>zh:</b> θ – 插管角度 (°) <b>pl:</b> θ – kąt rurki (°) <b>cs:</b> θ – Úhel trubice (°) <b>sk:</b> θ – Úhol kanyly (°) <b>hu:</b> C – tubus szöge (°) <b>el:</b> θ – Γωνία σωλήνα (°) <b>tr:</b> θ – Tüp Açısı (°) <b>th:</b> θ – ทüp Ang (°) <b>ko:</b> θ – 튜브 각도(°) <b>id:</b> θ – Sudut Tabung (°) <b>ar:</b> θ – زاوية الأنبوب (°)</p>
	<p><b>en:</b> Disposable Inner Cannula Length (mm) <b>fr:</b> Longueur de la canule interne jetable (mm) <b>de:</b> Länge der Innenkanüle zum Einmalgebrauch (mm) <b>nl:</b> Lengte wegwerpbare binnencanule (mm) <b>it:</b> Lunghezza controcanula monouso (mm) <b>es:</b> Longitud de la cánula interna desechable (mm) <b>sv:</b> Längd på innerkanyl för engångsbruk (mm) <b>da:</b> Længden på inderkanylen til engangsbrug (mm) <b>no:</b> Lengde på gjenbrukbar innerkanyl (mm) <b>fi:</b> Kertakäyttöisen sisäkanylin pituus (mm) <b>pt:</b> Comprimento da cânula interna descartável (mm) <b>ru:</b> Длина одноразовой внутренней канюли (мм) <b>zh:</b> 一次性内套管长度 (mm) <b>pl:</b> Długość jednorazowej kanuili wewnętrzej (mm) <b>cs:</b> Délka jednorázové vnitřní kanyly (mm) <b>sk:</b> Dĺžka jednorazovej vnútornej kanyly (mm) <b>hu:</b> Eldobható belső kanül hossza (mm) <b>el:</b> Μήκος αναλώσιμης εσωτερικής κάνουλας (mm) <b>tr:</b> Tek Kullanımlık İç Kanül Uzunluğu (mm) <b>th:</b> ความยาวของท่อใส่ยาแบบใช้ครั้งเดียว (ม.) <b>ko:</b> 일회용 내부 캐눌라 길이(mm) <b>id:</b> Panjang Kanula Bagian Dalam Sekali Pakai (mm) <b>ar:</b> طول القنية الداخلية المخصصة للاستخدام مرة واحدة (مم)</p>

	<p><b>en:</b> Test Inflation Air Volume (cc); syringe not included <b>fr:</b> Volume d'air de gonflage de test (cc) ; seringue non fournie <b>de:</b> Befüllungstest-Luftvolumen (ml); Spritze nicht inbegriffen <b>nl:</b> Luchtinflatievolume voor test (ml); spuit niet inbegrepen <b>it:</b> Volume aria di gonfiaggio per il test (cc); siringa non inclusa <b>es:</b> Prueba de volumen de aire de inflado (cc); jeringa no incluida <b>sv:</b> Luftvolym för testfyllning (cc), spruta ingår inte <b>da:</b> Testluftoppumpningsvolumen (cc), sprøjte medfølger ikke <b>no:</b> Luftvolum (cc) for testoppblasing, sprøyte ikke inkludert <b>fi:</b> Testitäytön ilmamäärä (cc), ruiskua ei toimiteta mukana <b>pt:</b> Volume do teste de insuflação de ar (cc); seringa não incluída <b>ru:</b> Тестовый объем накачивания воздухом (куб. см); шприц не входит в комплект <b>zh:</b> 检查用充气体积 (cc) ; 不包括注射器 <b>pl:</b> Objętość napełnienia testowego powietrzem (ml); brak strzykawki w zestawie <b>cs:</b> Zkušební plnicí objem vzduchu (cc); bez injekční stříkačky <b>sk:</b> Skúšobný objem nafuknutia vzduchom (cc), striekačka nie je súčasťou balenia <b>hu:</b> Teszteljesítmény levegőtérfogat (cm<sup>3</sup>); a fecskendő nincs mellékelve <b>el:</b> Δοκιμαστικός όγκος αέρα διαστολής (cc), δεν περιλαμβάνεται σύριγγα <b>tr:</b> Şişirme Havası Test Hacmi (cc); şiริงga dahil değildir <b>th:</b> ปริมาตรอากาศสำหรับทดสอบ (cc) <b>ko:</b> 시험 팽창 공기 볼륨 (cc), 주사기 불포함 <b>id:</b> Volume Udara Pengembungan Uji (cc); tidak termasuk alat suntik <b>ar:</b> حجم حواء الفتح (cc) ; المحقن غير مشمول</p>
	<p><b>en:</b> Cuff Resting Diameter (mm) <b>fr:</b> Diamètre du ballonnet au repos (mm) <b>de:</b> Cuff-Durchmesser (in Ruhe) (mm) <b>nl:</b> Diameter niet-opgeblazen cuff (mm) <b>it:</b> Diametro a riposo della cuffia (mm) <b>es:</b> Diámetro del manguito en reposo (mm) <b>sv:</b> Kuffens diameter (mm) <b>pt:</b> Diâmetro do manguito em repouso (mm) <b>ru:</b> Манжетанягсдiameter (mm) <b>no:</b> Mansjettens hvilediameter (mm) <b>fi:</b> Kuffin läpimitta tyhjänä (mm) <b>pl:</b> Diámetro do manguito em repouso (mm) <b>da:</b> Manøjet diameter i ro (mm) <b>zh:</b> 套囊静息直径 (mm) <b>pt:</b> Médida do manguito em repouso (mm) <b>cs:</b> Průměr manžety v klidu (mm) <b>sk:</b> Priemer manžety v klude (mm) <b>hu:</b> A mandzsetta felújítatlan állapotbeli átmérője (mm) <b>el:</b> Διάμετρος αερόβαλαίου πριν το φορέσιμό (mm) <b>tr:</b> Şişmemiş Kaf Çapı (mm) <b>th:</b> เส้นผ่าศูนย์กลางของสายรัดข้อเท้าก่อนใช้ (mm) <b>ko:</b> 커프 거치대 직경 (mm) <b>id:</b> Diameter Istirahat Manset (mm) <b>ar:</b> القطر الثابت الكفة (مم)</p>
	<p><b>en:</b> Loading Dilator <b>fr:</b> Dilatateur de chargement <b>de:</b> Ladedilatator <b>nl:</b> Laaddilatator <b>it:</b> Dilatatore di inserimento <b>es:</b> Dilataador de carga <b>sv:</b> Dilator för införande <b>da:</b> Indføringsdilatator <b>no:</b> Innføringsdilatator <b>fi:</b> Asetusdilatattori <b>pt:</b> Dilataador de carregamento <b>ru:</b> Установочный дилататор <b>zh:</b> 负载扩张器 <b>pl:</b> Rozszerzacz do zakładania <b>cs:</b> Zaváděcí dilatátor <b>sk:</b> Zavádzací dilatátor <b>hu:</b> Behelyezési dilatátor <b>el:</b> Διαστολέας τοποθέτησης <b>tr:</b> Yüklem Dilatörü <b>th:</b> เครื่องขยายรูปร่างของสายรัดข้อเท้า <b>ko:</b> 부하 확장기 <b>id:</b> Dilator Pemuatan <b>ar:</b> موضح التحميل</p>
	<p><b>en:</b> Obturator <b>fr:</b> Introducteur <b>de:</b> Obturator <b>nl:</b> Obturator <b>it:</b> Otturatore <b>es:</b> Obturador <b>sv:</b> Obturator <b>da:</b> Obturator <b>no:</b> Obturator <b>fi:</b> Obturaattori <b>pt:</b> Obturador <b>ru:</b> Обтуратор <b>zh:</b> 插管芯 <b>pl:</b> Obturator <b>cs:</b> Obturátor <b>sk:</b> Obturátor <b>hu:</b> Obturátor <b>el:</b> Εμπυροεικό <b>tr:</b> Obtüratör <b>th:</b> ใส่วัดรูปร่าง <b>ko:</b> 관개 <b>id:</b> Obturator <b>ar:</b> الفتحة</p>
	<p><b>en:</b> Neck Strap <b>fr:</b> Attache cervicale <b>de:</b> Halteband <b>nl:</b> Nekband <b>it:</b> Fettuuccia di fissaggio <b>es:</b> Cinta para el cuello <b>sv:</b> Halsrem <b>da:</b> Halsrem <b>no:</b> Halsstropp <b>fi:</b> Kiinnitysnauha <b>pt:</b> Correia para pescoço <b>ru:</b> Шейная лента <b>zh:</b> 颈带 <b>pl:</b> Pasek szyjny <b>cs:</b> Krční popruh <b>sk:</b> Popruh na krk <b>hu:</b> Nyakszij <b>el:</b> Ταβία στρίψινης στο λαιμό <b>tr:</b> Boyun Kayışı <b>th:</b> สายรัดคอ <b>ko:</b> 목 끈 <b>id:</b> Pengikat Leher <b>ar:</b> سربط الرقبة</p>

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**Shiley™****Adult Flexible Evac Tracheostomy Tube with TaperGuard™ Cuff  
Disposable Inner Cannula****Instructions for Use**

---

**Description**

The Shiley™ Adult Flexible Evac Tracheostomy Tube with TaperGuard™ Cuff Disposable Inner Cannula is a single patient use dual cannula tracheostomy tube. It contains a 15 mm connector on the outer cannula for use with standard ventilation and anesthesia equipment.

The tracheostomy tube has a clear outer cannula constructed of thermosensitive biocompatible polyvinyl chloride with a radiopaque line. A pliable neck plate allows conformity to individual neck anatomies and contains two suture holes. The tapered tracheostomy tube tip is compatible with percutaneous insertion. The smooth, hollow, and tapered tip obturator facilitates insertion and can be used with or without a guide wire. One disposable inner cannula is included in the package. The disposable inner cannula with embossed size is translucent for easy inspection and allows for airway maintenance. This tracheostomy tube can be used with or without the disposable inner cannula, however, the use of the inner cannula is recommended.

The tracheostomy tube has a taper-shaped, low-pressure cuff to ease insertion and provide an effective air and fluid seal. The cuff inflation line has a pilot balloon with a luer check valve for inflating and deflating the cuff, and monitoring the cuff pressure.

A clear suction line is attached above the neck plate that terminates to a lateral port above the cuff to permit the suctioning of subglottic secretions. The evacuation port is as close to the cuff as possible to avoid residual secretions and mucosal wall contact. The distal end of the clear suction line contains a yellow connector with a tethered cap for use with suction devices. The suction line is integrated into the wall of the outer cannula.

The tracheostomy tube and accessories with the exception of the cotton tape tie are composed of plastics which have not been made with natural rubber latex material or DEHP. They are safe to use for X-ray, ultrasound, Positron Emission Tomography (PET), and Radiation Therapy. The tracheostomy tube and accessories have been sterilized with ethylene oxide.

The sizing dimensions are listed on page 2 and refer to: the inside diameter (ID) of the tube or inner cannula at its narrowest point; the outside diameter (OD) of the outer cannula; and the length measured as the distance from the neck plate to the distal tip on the tube or inner cannula center line.

**Contents of Package:** 1 tracheostomy tube with cuff and pliable neck plate; 1 disposable inner cannula; 1 obturator; 1 neck strap. Reference page 2 for sizes and dimensions.

**Accessories Available Separately:** Sizes 6.5, 7.0, 7.5, 8.0, and 8.5 are compatible with the Percutaneous Dilatational Loading Dilator from the COOK Critical Care Ciaglia Percutaneous Introducer Trays. See chart on page 2 for tracheostomy and dilator sizing compatibility.

**Indications for Use**

The Shiley™ Adult Flexible Evac Tracheostomy Tube with TaperGuard™ Cuff is intended to provide tracheal access for airway management and subglottic secretion management in adult patients. It is also intended for use with percutaneous dilatational tracheotomy (PDT) procedures.

**Contraindications**

There are no known contraindications associated with this product.

**WARNINGS**

- **Avoid contact of tube with laser surgery beams or electrosurgery electrodes in the presence of oxygen enriched gas or nitrous oxide because this could result in combustion of the tube with harmful burns and emission of corrosive and toxic combustion products, including hydrochloric acid (HCl).**
- **This device must only be used by clinicians trained in its use.**

- Sterile only if the protective tray and cover are not open, damaged, or broken. Do not resterilize.
- Contents of this package cannot be resterilized by the user and are intended for single patient use only. Attempts to resterilize this device may result in product failure.
- Prior to use, handle with due care to prevent user contamination.
- Use expert clinical judgement when selecting a tracheostomy tube for a patient. The inner diameter, outer diameter, and length must be considered. If the inner diameter is too small, it will increase the resistance through the tube and make airway clearance more difficult. An outer diameter too small will increase the cuff pressure required to seal the trachea or lead to persistent air leak. An outer diameter too large will be difficult to pass through the stoma. A tube too short can obstruct the tube's distal end against the posterior tracheal wall or cause tracheal stenosis.
- Tube obstruction, false passage, or persistent air leak is more common when tube is used in patients with increased pretracheal tissue or fat.
- Test each tube's cuff, pilot balloon, and valve by inflation prior to use. If dysfunction is detected in any part of the inflation system, the tube should not be used and should be returned to the manufacturer for investigation.
- Inflation volumes listed on page 2 are for cuff test purposes only.
- Inflation of the cuff by "feel" alone or by using a measured amount of air is not recommended. The pilot balloon is only intended to indicate the presence of pressure or vacuum in the cuff and is not intended to provide an indication of pressure level.
- Cuff pressure should be monitored and adjusted routinely and should not exceed 30 cm H<sub>2</sub>O (3 kPa). Over-inflation of the cuff may cause tracheal damage and inhibit ventilation. Under-inflation with cuff pressures less than 20 cm H<sub>2</sub>O (2 kPa) may result in aspiration of subglottic secretions.
- Diffusion of nitrous oxide, oxygen, or air can alter cuff volume and pressure.
- Altitude, temperature, body position, esophageal medical instruments, and abdominal insufflation can alter cuff pressure.
- Suctioning and tube care should be performed frequently enough to prevent the buildup of secretions.
- Use of proper humidification of inspired air, use of an inner cannula, and attentive nursing care will minimize potential secretion build up within the tube.
- When using an inner cannula, it can only be replaced by a Shiley™ inner cannula of the same size. It is recommended that a spare disposable inner cannula be with the patient at all times.
- A replacement tracheostomy tube and obturator should be kept at bedside.
- Do not apply greater than 150 mmHg (20 kPa) of intermittent suction pressure to the subglottic port to avoid tracheal injury.
- Ventilator and oxygen-delivery tubing should be attached without exerting traction on the tracheostomy tube because repeated motion of the tubing can lead to skin breakdown, stenosis and / or tracheal tissue damage.
- Cap the suction line when not attached to the suction device to avoid the risk of infection.
- A fully seated obturator must always be used during tube changes or reinsertion during unplanned decannulations to avoid tracheal damage. Always keep replacement device and obturator at an accessible location near the patient.
- Prior to the removal of the tube from the patient, disconnect suction source, and completely remove all air from the cuff to prevent damage to the trachea and stoma.
- When using a water soluble lubricant, verify that the lubricant does not enter and occlude the tube lumen which may restrict ventilation.
- To prevent patient breathing difficulty and/or suffocation, the tracheostomy tube cuff must be completely deflated before placing a Shiley™ speaking valve.

## CAUTIONS

- Federal law (USA) restricts this device to sale by or on the order of a physician.

- Duration of patient use of the tracheostomy tube and obturator should not exceed twenty-nine (29) days because the manufacturer has not substantiated the use of these devices beyond 29-days. Decisions about tracheostomy tube changes should be made by the responsible physician or designate using current medical guidelines and judgment.
- Visually inspect device and accessories for signs of damage. If signs of damage, discard and replace.
- Prior to use, check that the obturator can be inserted and removed without difficulty.
- Clinicians placing this device percutaneously should be competent in the procedure prior to using this tube. To ease placement of the tracheostomy tube, generously lubricate the surface of the loading dilator.
- When deflating the cuff, position it away from the distal tip toward the neck plate. When tapering the cuff, do not use sharp instruments that would damage the cuff, such as forceps or hemostats.
- Avoid pulling or manipulation of the inflation and suction lines in order to maintain effective cuff and evac performance.
- Prevent lint or other particulates from entering the luer valve of the pilot balloon.
- Extremes of tube movement may result in mucosal abrasion and transient disruption of cuff seal.
- The use of Lidocaine aerosol has been associated with the formation of pin holes in PVC cuffs. Lidocaine hydrochloride solution has been reported not to have this effect.
- Do not clean or reuse device or accessories..
- A lint-free gauze pad soaked with Sterile Saline can be used to wipe the device; do not use solutions or chemical agents.

### **Use During MRI**

Non - clinical testing has demonstrated that the device is MR Conditional. This device includes a luer valve that contains a metal spring that is external to the patient. A patient with this device can be safely scanned in a MR system meeting the following conditions:

- Static magnetic field of 1.5 T or 3.0T, maximum spatial field gradient of 1,900 gauss/cm (19T/m) maximum MR system reported whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode).
- Under the scan conditions defined above, the metal spring in the luer valve is expected to produce a maximum temperature rise of 0.2+–1.5C after 15 minutes of continuous scanning.
- MR image quality may be compromised if the imaging area of interest is close to the position of the luer valve and may require optimization of MR imaging parameters.
- In non-clinical testing using a 1.5T MRI System, the image artifact caused by the luer valve extends approximately 72 x 44 mm with a spin echo sequence and 103 x 43 mm with a gradient echo sequence. In non-clinical testing using a 3T MRI System, the image artifact caused by the luer valve extends approximately 30 x 27 mm with a spin echo sequence and 118 x 43 mm with a gradient echo sequence.
- During an MRI scan the luer valve should be secured away from the area of interest to prevent movement and image distortion.

### **Tube Preparation and Cuff Inflation Test**

1. Per the physician's order, select the appropriate tracheostomy tube size.
2. Use sterile technique for the handling and placement of the tracheostomy tube.
3. Test the cuff and inflation system for leakage before inserting the tube, as follows:
  - a. Inflate the cuff with the volume of air indicated on page 2.
  - b. Immerse the tube in sterile saline and observe for air leakage. If any leak is noted, do not use the tube and return to the manufacturer for investigation.
  - c. When deflating the cuff, position cuff away from the distal tip toward the neck plate. This will ease insertion and guard against cuff perforation from sharp edges of cartilage.

### **Tube Insertion Using the Obturator or a Percutaneous Loading Dilator**

1. Select the compatible loading dilator for the tube size listed on page 2 for percutaneous dilatational tracheotomy (PDT) procedures. (The loading dilator is available separately.)
2. Insert the obturator or compatible loading dilator into the tracheostomy tube.

- If using the **obturator**, ensure it is fully seated. Apply a thin layer of water soluble lubricant to the outside of the tracheostomy tube, cuff, and protruding portion of the obturator. If using a guide wire, after insertion of the wire into the patient's stoma thread the wire all the way through the obturator starting at the distal tip.
  - If using a **loading dilator for a PDT procedure**, generously lubricate the surface of the appropriately sized loading dilator and insert into the tracheostomy tube. Ensure that the tapered section of the loading dilator clears the distal tip of the tracheostomy tube by approximately 2 cm and that the transition of the dilator with the tip of the tube is smooth.
3. Insert the tube into the patient's trachea and immediately remove the loading obturator/dilator and guide wire. The replacement device and obturator should be stored in an accessible location near the patient if needed for tracheostomy tube re-insertion.
  4. Inflate the cuff to an intracuff pressure of 25–30 cm H<sub>2</sub>O (2.5-3 kPa).
  5. Secure the tracheostomy tube to the patient using the cotton neck strap provided such that only one finger may be inserted between the tie and the skin of the neck before dressings are applied.
  6. Insert the Disposable inner cannula into the tracheostomy tube by squeezing the proximal end of the inner cannula and lightly pressing the end of the inner cannula to lock into position.

### Subglottic Suction

1. Check that the cuff is inflated to an intracuff pressure of 25–30 cm H<sub>2</sub>O (2.5-3 kPa), then attach the suction line to a suction device using intermittent suction of 100-150 mmHg (13.3-20 kPa). Do not exceed 150 mmHg (20 kPa).
2. Monitor intracuff pressure every 4 hours and maintain at 25–30 cm H<sub>2</sub>O (2.5-3 kPa).
3. Every 4 hours check that the suction lumen is patent by disconnecting the tube suction line and injecting 2 mL of air or 5 mL of sterile 0.9% saline solution with a 10 mL syringe. Reconnect the suction line to the suction device. If no saline or secretions are recovered or if saline leaks from the stoma, discontinue suction and notify the attending physician for evaluation of proper tube size or placement.
4. Monitor the suction pressure at least every 4 hours. If the suction pressure increases above the set pressure, check the system as in the previous steps.
5. Cap the yellow connector when not attached to the suction device to prevent contamination and possible infection of the tube.

### Changing the Disposable Inner Cannula

1. A color-coding system is used to appropriately match replacement DIC with the corresponding same-sized tube. This color-coding can be found on the outer box label on the disposable inner cannulas and on the flange of the tube.
2. While holding on to the tracheostomy tube's 15 mm connector, gently twist and pull the breathing circuit connector to disconnect. Remove the disposable inner cannula by squeezing the end of the inner cannula and slowly withdrawing it.
3. Insert a new disposable inner cannula by squeezing the proximal end of the inner cannula and lightly pressing into locked position.
4. Reconnect the breathing circuit.

### Tube Removal

1. Remove the neck strap from around the patient's neck.
2. Deflate the cuff using a syringe to slowly withdraw air from the luer valve of the inflation line. Observe that the pilot balloon is completely collapsed. Deflating the cuff ensures that the cuff passes through the stoma with minimum resistance.
3. Slowly remove the tracheostomy tube from the patient.
4. Dispose of the tracheostomy tube and accessories in accordance with current medical standards and applicable national regulations for biologically hazardous waste.

### Adverse Events

The following adverse events have been reported to be associated with the use of tracheostomy tubes during the tracheotomy procedure, during cannulation, or after decannulation, listed in order of frequency and severity: Infection / aspiration, tissue trauma, dislodgement, minor bleeding, false tract, partial/full obstruction, severe tissue trauma / hemorrhage, change in vital signs / arrest, gas trapping, dysphagia, increased work of breathing, skin/stoma breakdown, vocal cord paralysis, difficult/multiple intubation, hypoxia, hypercapnia, difficult weaning, airway fire, foreign body, toxicity, barotrauma, image distortion if used in MRI scanner.

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Not made  
with DEHP

**DEHP**

Do not use if  
package is opened  
or damaged

MR Conditional

0123

Part No. PT00058806 Rev C 2019-05

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. <sup>TM</sup>\* brands are trademarks of their respective owners. Other brands are trademarks of a Covidien company.

© 2017 Covidien.

 Covidien llc, 15 Hampshire Street,  
Mansfield, MA 02048 USA.

**EC REP** Covidien Ireland Limited,  
IDA Business & Technology Park,  
Tullamore, Ireland.

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