Shiley™

Tracheostomy Tube

REF 50XLTC D

REF 50XLTCP

REF 50XLTIN

REF 50XLTUD

REF 50XLTUP

REF 60XLTC D

REF 60XLTCP

REF 60XLTIN

REF 60XLTUD

REF 60XLTUP

REF 70XLTC D

REF 70XLTCP

REF 70XLTIN

REF 70XLTUD

REF 70XLTUP

REF 80XLTC D

REF 80XLTCP

REF 80XLTIN

REF 80XLTUD

REF 80XLTUP

en

Instructions For Use

fr

Mode d'emploi

de

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it

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es

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pt

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nl

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sv

Bruksanvisning

da

Brugsanvisning

no

Bruksanvisning

fi

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el

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XX = First 2 Digits of Product Code Designated to the Inner Diameter.

Table 2


REF								
				A (mm)	B (mm)	C (mm)		
50XLTCD	5.0	9.6	90	48.0	5.0	37.0	29.2	25
50XLITUD	5.0	9.6	90	48.0	5.0	37.0		
50XLITCP	5.0	9.6	90	33.0	20.0	37.0	29.2	25
50XLITUP	5.0	9.6	90	33.0	20.0	37.0		
60XLTCD	6.0	11.0	95	49.0	8.0	38.0	31.3	26
60XLITUD	6.0	11.0	95	49.0	8.0	38.0		
60XLITCP	6.0	11.0	95	34.0	23.0	38.0	31.3	26
60XLITUP	6.0	11.0	95	34.0	23.0	38.0		
70XLTCD	7.0	12.3	100	49.0	12.0	39.0	35.0	29
70XLITUD	7.0	12.3	100	49.0	12.0	39.0		
70XLITCP	7.0	12.3	100	34.0	27.0	39.0	35.0	29
70XLITUP	7.0	12.3	100	34.0	27.0	39.0		
80XLTCD	8.0	13.3	105	50.0	15.0	40.0	35.0	29
80XLITUD	8.0	13.3	105	50.0	15.0	40.0		
80XLITCP	8.0	13.3	105	35.0	30.0	40.0	35.0	29
80XLITUP	8.0	13.3	105	35.0	30.0	40.0		

Table 3

REF Cat. No.	REF Cat. No.		
XLT Disposable Inner Cannula	For use with XLT Tracheostomy Tube	Inner Cannula I.D. (mm)	Inner Cannula Length (mm)
50XLTIN	50XLTUP 50XLTCP 50XLTUD 50XLTCB	5.0	90
60XLTIN	60XLTUP 60XLTCP 60XLTUD 60XLTCB	6.0	95
70XLTIN	70XLTUP 70XLTCP 70XLTUD 70XLTCB	7.0	100
80XLTIN	80XLTUP 80XLTCP 80XLTUD 80XLTCB	8.0	105

Shiley™ Tracheostomy Tube

This instruction insert applies to the following Shiley products:

Tracheostomy Tube XLT Cuffed Distal with Disposable Inner Cannula
 Tracheostomy Tube XLT Cuffed Proximal with Disposable Inner Cannula
 Tracheostomy Tube XLT Cuffless Distal with Disposable Inner Cannula
 Tracheostomy Tube XLT Cuffless Proximal with Disposable Inner Cannula
 Inner Cannula XLT Disposable

Identification of a substance that is contained or present within the product or packaging.

Identification of a substance that is not contained or present within the product or packaging.

Cuff resting diameter

This product contains DEHP. When used as indicated, very limited exposure to trace amounts of DEHP may occur. There is no clear clinical evidence that this degree of exposure increases clinical risk. However, in order to minimize risk of DEHP exposure in children and nursing or pregnant women, this product should only be used as directed.

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

Directions for Use

Indications

The devices are intended for use in providing tracheal access for airway management.

WARNINGS

- **The Shiley Tracheostomy Tube XLT and obturator are single patient use medical devices and usage should not exceed twenty-nine (29) days. Frequent and routine changes of the tracheostomy tube and accessories are recommended and should be evaluated by the qualified physician.**
- **The Shiley Inner Cannula XLT, neck strap and Shiley TTH Tracheostomy Tube Holder are disposable medical devices and are designed for single use and should not be cleaned or reused.**
- **Carefully read these instructions prior to use. These instructions are applicable only to the Shiley Tracheostomy Tube XLT, Shiley Inner Cannula XLT.**
- **The products packaged in the tray inside the box are sterile unless the package is opened or damaged. Do not resterilize. Do not expose to temperatures above 120°F (49°C).**
- **The Shiley Inner Cannula XLT provided in this package can only be replaced by a Shiley Inner Cannula XLT of the same size and catalog number.**
- **A Shiley Inner Cannula XLT must be in place to mechanically ventilate the patient.**

- **During and after attachment of respiratory or anesthesia tubing connectors to the Shiley™ Inner Cannula XLT, avoid application of excessive rotational, linear, or rocking forces on the tubing and/or connectors to prevent damage to the tracheostomy tube.**
- **Because of its special design, the Shiley Inner Cannula XLT requires careful handling during insertion and removal. When inserting an inner cannula, use a smooth, uniform and consistent motion, using the least amount of force necessary for complete insertion.**
- **If you are unable to insert or remove the inner cannula, contact a physician immediately.**
- **Under no circumstances should more than 25 mm of mercury (Hg) air pressure be used to inflate the cuff. Over-inflation of the cuff may cause tracheal damage and may inhibit ventilation.**
- **As a further precaution for ventilator-dependent patients, cuff inflation should be checked on a regular basis and replacement tracheostomy tubes and inner cannula should be kept at the patient's bedside**
- **Do not use the Shiley Tracheostomy Tube XLT and Shiley Inner Cannula XLT in procedures which will involve the use of LASER or an electro-surgical active electrode in the immediate area of the devices. Contact of the beam or electrode with the tracheostomy tube and inner cannula, especially in the presence of oxygen-enriched or nitrous oxide containing mixtures could result in the rapid combustion of the tube with harmful thermal effects and toxic emissions.**
- **During an MRI scan the pilot balloon should be secured near the Y connector of the ventilator circuit at least 3 cm from the area of interest to prevent movement and image distortion.**

CAUTIONS

- Precautions should be taken when discarding any of the tracheostomy tube devices and accessories. Disposal of medical devices should be made in accordance with applicable national regulations for biologically hazardous waste.
- The Shiley Inner Cannula XLT, neck strap and Shiley TTH Tracheostomy Tube Holder are designed for single use and should not be cleaned or reused.
- Patients in the home care environment should be carefully instructed by the Home Health Care Provider in the proper use and handling of the tracheostomy tubes and accessories.
- To ease insertion and to guard against cuff perforation from sharp edges of cartilage, the cuff should be tapered back. This can be accomplished by first inflating the cuff. Then gently move the cuff away from the distal tip of the outer cannula towards the swivel neck plate as the residual air is removed by deflation. Do not use any sharp instruments such as forceps or hemostats that would damage the cuff when tapering it.
- When using a water soluble lubricant, verify that the lubricant does not enter and occlude the tube lumen, which may restrict ventilation. Petroleum-based lubricants should not be used.
- The obturator should be cleaned (Refer to Table 4) and then stored in an accessible location near the patient should it be needed.
- To avoid applying pressure against the patient, the neck plate may be stabilized with your free hand during the locking procedure.
- Verify that the twist-lock connector engages securely. If parts become worn or loose, immediately report this to the physician for prompt replacement of the tracheostomy tube.
- This product is composed of soft materials to conform to tracheal tissue for performance and patient comfort. Simple precautions in handling of the Shiley

Tracheostomy Tube XLT with Cuff during insertion and while in place will facilitate proper function and minimize tears and breaks in the cuff. Avoid pulling or manipulation of the inflation line, as it is designed to conduct and hold air as part of the cuff inflation system. It is recommended that the inflation line be maintained in a position allowing for patient mobility, without placing tension on the line-to-outer cannula junction. Prevent lint or other particulates from entering the luer valve of the pilot balloon.

- Use the TTH Tracheostomy Tube Holder in supervised care settings only. Pediatrics, confused patients and other individuals who require special care should be monitored to prevent unintentional detachment.

NOTE: Shiley™ Tracheostomy Tube XLT, Shiley Inner Cannula XLT, Shiley TTH Tracheostomy Tube Holder and Accessory Products are composed of medical grade materials which do not contain latex.

Contents of Packages - Refer to Table 1.

Description - Refer to Table 2 and 3 for Sizing Dimensions.

NOTE: The sizing dimensions refer to the I.D. as the inside diameter of the inner cannula at the smallest location. O.D. is the outside diameter of the outer cannula (not including cuff features). Length is the distance from the patient's side of the neck plate to the distal tip on the tube's centerline.

The Shiley Tracheostomy Tubes XLT are biocompatible, double cannula tracheostomy tubes with disposable Inner Cannulae. The outer cannula is constructed of thermosensitive polyvinyl chloride with a radiopaque line to assist in radiographic visualization. A soft swivel neck plate allows conformity to individual neck anatomies. The Shiley Inner Cannula XLT with 15 mm connector, is translucent for easy inspection and can be used for mechanical ventilation with standard respiratory equipment. The smooth, rounded tip obturator facilitates insertion. The Shiley Tracheostomy Tubes XLT are available in four sizes: 5.0mm, 6.0mm, 7.0mm, and 8.0mm I.D. Additional Shiley Inner Cannula XLT in the same sizes are available separately in a 10 pack carton.

The cuffed Shiley Tracheostomy Tube XLT has a thin wall, high volume, low pressure cuff to minimize pressure in the trachea. When inflated, the cuff conforms to the natural shape of the trachea, providing a seal at low intracuff pressure. The cuff inflation line has a luer valve with an integral pilot balloon to indicate cuff inflation.

The Shiley TTH Tracheostomy Tube Holder is a comfortable two piece, one size fits all, holder used to secure all tracheostomy tubes. Made of soft breathable cotton lined foam laminated with Velcro® fasteners.

Velcro® is a registered trademark of Velcro Inc., USA

Instructions for use

Tube Preparation

1. Selection of tracheostomy tube size is left to the discretion of the physician. Patients in the home care environment should be carefully instructed by a Home Health Care Provider in the proper use and handling of this device.

Pre-Insertion Cuff and Inflation Test

NOTE: Refer to (Table 2) for Leak Test Inflation Volumes. Inflation volumes are for test purposes only. Consult the physician or Home Health Care Provider for the appropriate inflation volume/pressure when the tube is positioned within the trachea.

2. With the cuffed Shiley™ Tracheostomy Tube XLT, the cuff and inflation system should be tested for leakage before inserting the tube. This test can be performed as follows: Inflate the cuff with the volume of air indicated in (Table 2). Then either observe for deflation over several minutes or immerse the tube in sterile saline and observe for air leakage. Completely deflate the cuff prior to insertion.

Insertion

3. Insert the obturator into the tracheostomy tube. The obturator should be fully seated before the tracheostomy tube is inserted into the patient. A thin film of water soluble lubricant can be applied to both the outer cannula, cuff, and protruding portion of the obturator to facilitate insertion.
4. After the tracheostomy procedure is performed, insert the tube into the patient's trachea. After the tube is properly in place, remove the obturator immediately.
5. The Shiley Inner Cannula XLT is thin walled to provide maximum airway diameter and has a special flexible design to allow it to assume the shape of the Shiley Tracheostomy Tube XLT in-situ.
 - (a) Always inspect the Shiley Inner Cannula XLT prior to use. If the inner cannula is damaged, creased or the lumen is compromised, discard and replace with a new inner cannula. Do not insert a damaged inner cannula.
 - (b) Insert the Shiley Inner Cannula XLT into the Shiley Tracheostomy Tube XLT and lock into position. The Shiley Inner Cannula XLT may be moistened with sterile saline or water soluble lubricant to facilitate insertion. To lock the Shiley Inner Cannula XLT in place, twist the locking ring clockwise until it stops. The Shiley Tracheostomy Tube XLT should always be used with the Shiley Inner Cannula XLT in place.

NOTE: The inner cannula is locked only when the indicator points of the inner cannula locking ring and tube outer cannula are aligned.

6. The Shiley Inner Cannula XLT can be replaced to maintain a patent airway without requiring the entire tracheostomy tube to be changed.
7. To remove the Shiley Inner Cannula XLT from the Shiley Tracheostomy Tube XLT, twist the locking ring counter-clockwise until it stops. Withdraw slowly and discard.

Cuff inflation following insertion


8. Inflate the low pressure cuff by injecting air into the luer valve of the inflation line. Selection of a cuff inflation and deflation procedure is left to the discretion of the clinician.
9. Secure the Shiley Tracheostomy Tube XLT to the patient using either the TTH Tracheostomy Tube Holder or neck strap provided.

Diagram 1

Diagram 2

Diagram 3

Diagram 4

- 1-2 **Thread narrow Velcro® fasteners through slots in flange and adhere to foam laminate.**
- 3-4 **Adjust bands to proper length around patient's neck, adhere wide Velcro® fastener (A) to band (B). Cut off excess.**

Cuff Deflation

10. Accumulated secretions above the cuff may need to be suctioned before deflating the cuff.

11. To deflate the low pressure cuff, withdraw air slowly from the luer valve of the inflation line using a syringe.

Removing Tube

12. Before removing the entire tracheostomy tube, fully deflate the cuff (if present) with a syringe. This will ensure that the cuff passes through the stoma with minimal resistance.

Cleaning

WARNINGS

- **Do not use solutions or chemical agents other than those recommended in Table 4 to clean any part of the Shiley™ Tracheostomy Tube XLT, as this may result in damage.**
- **The Shiley Inner Cannula XLT and Shiley TTH Tracheostomy Tube Holder should not be cleaned or reused.**

Table 4

PART	RECOMMENDED CLEANING
<ul style="list-style-type: none">• Disposable inner cannula• TTH Tracheostomy Tube Holder• Neck strap	DO NOT CLEAN OR REUSE.
<ul style="list-style-type: none">• Cuffed Outer Cannula	Rinse Gently with Sterile Saline. Air Dry.
<ul style="list-style-type: none">• Cuffless Outer Cannula• Neck Plate• Obturator	<ol style="list-style-type: none">1. Clean with hydrogen peroxide (half-strength), sterile normal saline, or water and mild detergent.2. After cleaning thoroughly rinse with sterile saline to remove any cleaning solution residues.3. Air dry.

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






Rx
ONLY

Do not
resterilize

Do not use if
package is opened
or damaged

Caution, consult
accompanying
documents

Upper temperature
limit

Part No. 10093381 Rev B 06/2013

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.

Velcro is a registered trademark of Velcro Industries B.V.

© 2011 Covidien.

 Covidien Inc, 15 Hampshire Street,
Mansfield, MA 02048 USA.

 Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

0123

www.covidien.com