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Intended use/Indication

The artificial nose (heat and moisture exchanger (HME)) is used for air conditioning the respiratory gases of spontaneously breathing patients with an endotracheal tube or tracheostomy tube. The integrated oxygen port is intended for the administration of additional oxygen. This product is a medical device.

User group

The user must possess the necessary qualifications, e.g. medical expertise or be trained in this respect and/or possess experience with the device.


The user must read and understand the Instructions for Use before use in order to use the device as intended.

Patient group

The artificial nose is intended for adult and paediatric patients.

Product description

REF	Type	Safety valve
43.009.05.705	301	no
43.009.05.706	302	yes

Type 302 with safety valve (B) for forced inspiration/expiration

The safety valve on Type 302 allows high tidal volumes for short periods and enables expectoration of larger volumes of secretions.

Technical specifications: (DIN EN ISO 9360-1, DIN EN ISO 9360-2)

- Tidal volume:** 50 - 1000 ml
- Dead space:** 17 ml
- Moisture loss:**
 - VT = 250 ml: 9.0 mg/l
 - VT = 500 ml: 13.2 mg/l
 - VT = 750 ml: 14.2 mg/l
 - VT = 1000 ml: 15.5 mg/l
- Pressure loss:**
 - Inspiration**
 - 0.5 l/s = 0.09 cm H₂O (~hPa)
 - 1.0 l/s = 0.26 cm H₂O (~hPa)
 - 1.5 l/s = 0.54 cm H₂O (~hPa)
 - Expiration**
 - 0.5 l/s = 0.11 cm H₂O (~hPa)
 - 1.0 l/s = 0.32 cm H₂O (~hPa)
 - 1.5 l/s = 0.34 cm H₂O (~hPa)
- Weight:** 5 grammes
- Length:** 29 mm

Connectors: 15 mm female for endotracheal tube and tracheostomy tube, 5 mm connector for oxygen connection tube

Contraindications

The artificial nose REF 43.009.05.705 should not be used in patients with excessive secretion from the airways. We are not aware of any other contraindications.

Material

Housing: polypropylene
HME medium: polyurethane
Safety valve: thermoplastic polymer
The device is latex-free.

Application

- Do not use if the packaging and/or the device are damaged or if a malfunction of the device is identified.
- The artificial nose is to be attached to the distal end of the endotracheal tube or tracheostomy tube.
- Check the artificial nose for a tight fit and leak tightness.
- If required, an oxygen connection tube can be connected to the oxygen port (A).
- The device must be checked regularly for fluid and secretion accumulation and changed immediately if necessary.
- The safety valve (B) of device variant Type 302 must always remain free and flexible
- This device must be changed every 24 hours.
- This device is intended for single use only. Re-use leads to multiplication and spreading of germs.

Notes

- Read the Instructions for Use carefully.
- The artificial nose should only be used for its intended purpose to ensure its safe use.
- Do not use the device in combination with active humidifiers or medication nebulisers. This can lead to a build-up of pressure and insufficient breathing of the patient.
- Do not use the device with other components such as filters or HMEs, as these increase the dead space and lead to an increase in compliance and resistance.
- If oxygen is supplied, the artificial nose is not suitable for smokers due to the risk of fire. Keep away from naked flames, sparks or other sources of ignition.
- The oxygen port (A) is only to be used for introducing oxygen.
- Sterilisation and cleaning (for re-use) can impair the function and mode of operation of the device and cause a hazard for the patient.
- The device must not be modified as this could lead to malfunctions and unforeseen risks such as injury to the patient and user or damage to property.

Disposal

The artificial nose must be disposed of in accordance with the in-house instructions and local regulations.

Storage/ambient conditions

Observe the permissible ambient temperature ranges when storing and using the artificial nose.

Additional information

All serious incidents occurring in connection with the medical device must be reported to the manufacturer and the responsible authorities.

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P. J. Dahlhausen & Co. GmbH